

## Set up Kohler illumination

Note: Only need to do this once at the beginning of imaging session!

IMPORTANT: Sample should be focused before hand

1. Close Field stop aperture (to left)
2. Open condenser aperture
3. Move condenser to height to superimpose octagon on sample
  - Can look through screen or eyepieces to do this
4. Look through eye piece to position octagon in center with condenser screws
  - Must be done through eye pieces
  - If light is too bright use filters at the top (Push in)

TIP: If you need to switch between eye and camera, use light path buttons

