



Fluorescence, Fluorescent microscopy and probes



Fluorescence in live cells








What is it?

Sir John Frederick William Herschel, 1854: Though perfectly transparent and colorless when held between the eye and the light, or a white object, it yet exhibits in certain aspects, and under certain incidences of the light, an extremely vivid and beautiful celestial blue colour, which, from the circumstances of its occurrence, would seem to originate in those strata which the light first penetrates the liquid.....


Fluorescence

Jablonski diagram






Polarization/Anisotropy



Why use Fluorescence probes?


- Sensitivity
- Specificity
- Analysis of location and quantity of a single component in a complex mixture
- Detection of small quantities of fluorophores and fluorescent objects below the resolution limit
- Environmental sensitivity
- Does not rely on physical properties of the specimen for contrast generation



Commonly Used Light Sources for Fluorescence Microscopy


- Mercury arc lamp
- Xenon arc lamp
- Metal Halide doped Hg-Xe Arc
- LASERs
(most often used for confocal, TIRF)

Light Sources: Mercury & Xenon Arc Lamps



Metal Halide doped Hg-Xe Arc

- Halogens decrease carbon deposits and slow deterioration of electrodes
- X-cite illuminator uses liquid light guide; minimal heat transfer



Laser lines

A line graph showing the relative intensity of laser lines for Argon, HeNe, and Krypton. The y-axis is relative intensity from 0.00 to 1.00. The x-axis is wavelength from 300 to 700 nm. Argon lines are shown in purple (around 400 nm), HeNe lines in green (around 546 nm), and Krypton lines in red (around 647 nm). A color bar at the bottom indicates the visible spectrum.

<http://www.repairfaq.org/sam/lasersam.htm>

A schematic diagram of a fluorescence microscope. It shows a 'Light Source' (1) emitting light through an 'Excitation filter' (2) and a 'Dichroic mirror' (2) to the sample. The sample emits light through the 'Dichroic mirror' (2) and an 'Emission filter' (2) to the eyepiece. A detailed inset shows the internal components: 'Exciter Filter', 'Dichroic Mirror', and 'Barrier Filter'.


4

Filters


- Need to reject excitation light completely
- Need to be transparent for emitted light
- Need to match spectra of dyes
- Spectra of dyes:

www.zeiss.com/micro
probes.invitrogen.com/resources/spectraviewer/
www.mcb.arizona.edu/ipc/fret/


Interference Filter Design (one cavity)



Interference Filter Design (multiple cavities)



Filter Terminology



Choose filters that maximize excitation and emission




Image from www.probes.com

Choose filters that maximize excitation and emission






Image from www.probes.com


Choose filters that separate fluorophores







Two different UV filter sets



Filter cube (after Ploem)








Types of fluorescent probes

- Immunofluorescence
- Fluorescent small molecules that bind specific cellular structures
 - DNA intercalating dyes (DAPI)
- Fluorescently labeled small molecules that bind specific cellular structures
 - Fluorescent phalloidin or taxol
- Fluorescently labeled proteins
- Fluorescent proteins (GFP)
- Genetically encoded tags binding fluorescent small molecules

The 'classic' dyes...



The Alexa Series Emission Spectra



Molecular Probes (www.probes.com)

Conjugation of organic dyes

Chemistry/Method

Amino groups (lysine): succinimidyl ester or isothiocyanate


Example:
Dynein driven gliding of microtubules
labelled with TMR on lysine residues.



Targets


- Antibodies: direct/indirect labeling (Label density)
- Proteins: labeling site unspecific
- Small molecules, i.e. phalloidin, taxol

Immunofluorescence




Quantum dots


- nanometre-scale crystals composed of atoms of an inorganic semiconductor material



Small dyes targeting specific cellular targets




DAPI
Hoechst 33258
Hoechst 33342
~20 fold enhancement



Mitotracker
Oxidized in mitochondria in
fluorescent compound



GFP-discovery




O. Shimomura
D.C. Prasher
M. Chalfie

<http://www.conncoll.edu/ccacad/zimmer/GFP-ww/GFP2.htm>

GFP structure





The 2004 palette of nonligomerizing fluorescent proteins



Roger Tsien


Other genetically encoded tags

Biarsenical-tetracysteine labeling




HaloTag

- Catalytically inactive mutant of a hydrolase that efficiently forms a covalent bond with HaloTag ligand (Promega)




Quenching and FRET





Use FRET as molecular ruler

$$E = R_0^6 / (R_0^6 + R^6)$$




Evidence for single dye pair FRET

When Cy5 bleaches, Cy3 emission recovers




FRAP



Ellenberg et al., 1997 (Lippincott-Schwartz lab)

(Fluorescence Recovery After Photobleaching)


FLIP



Lippincott-Schwartz lab


(Fluorescence Loss in Photobleaching)

Photo-activation (PA-GFP)




Excitation spectrum
before (filled), and
after (open)
photoactivation

Photo-activation (PA-GFP)




Koehler illumination



Lab samples:

- Drosophila S2 cells with various GFP-fusion proteins
- Stained tissue culture cells (multi-channel fluorescence)
- Fluorescent beads (visualize point spread function, registration shift)

THANKS!



Jennifer Waters

www.micro.magnet.fsu.edu

www.mcb.arizona.edu/ipc/fret/