


# Fluorescence Microscopy


Kurt Thorn  
NIC



# Why fluorescence?

- High contrast
  - Signal against dark background
- Highly specific, multi-color labeling
  - GFP etc.
  - Antibodies
- Live imaging
  - GFP etc.
- Quantitative
- Sensors for [Ca], pH, ...


# What is fluorescence?



Emission light is longer wavelength (lower energy) than excitation light


# Jablonski diagram

(Molecular energy diagram)




# Jablonski diagram

(Molecular energy diagram)




# Jablonski diagram

(Molecular energy diagram)




# Singlet and Triplet States



Spin flips are “dipole forbidden” → unlikely → long triplet lifetime



# Fluorescence Spectra



Alexa 488



# The Epifluorescence Microscope









Figure 2




Bandpass




Longpass






$OD = -\log_{10}(\%T)$



# Interference filters



Interference


→ Wavelength-  
Dependent  
transmission  
& reflection

→ filter

Filter makers :  
Chroma  
Semrock  
Omega


# Filter spectrum is angle dependent

Generally bad



# Filter spectrum is angle dependent



But can be good – Semrock VersaChrome



# Blocking

Interference filters have finite stop bands

Unblocked bandpass interference filter




Semrock 697/75

To block unwanted transmission from UV to IR, filter makers add *absorption glass* to the filter.

Often excitation filters are blocked, but emission filters *unblocked*.

→ Red autofluorescence or room light may get through your blue emission filter






## Filter cube nomenclature

- Chroma labels filters as center wavelength / passband (e.g. D350/50x)
- Dichroics are labeled by cut-on wavelength (e.g. 505DCLP)
- Nikon filters use a letter to specify illumination wavelength – e.g. UV, B, G, R
- Letters afterward specify emission profile – e.g. UV-2A vs UV-2E/C

# Exploding the filter cube

Emission  
Filter




Dichroic  
Mirror

Excitation Filter  
(off slide)

# Filter schemes

## Single wavelength sets


- Most efficient
- Best separation
- Very slow to change  $\lambda$



## Multi-band filters

### • Multi-band everything

- See all colors at once
- For color cameras
- Bad crosstalk




Chroma triple Pinkel set

### • “Pinkel” scheme

Multi-band dichroic  
Multi-band emitter  
Single- $\lambda$  exciters



- Excitation filter wheel
- Separate image at each wavelength
- Better separation



### • “Sedat” scheme

Multi-band dichroic  
single-band emitters  
Single- $\lambda$  exciters

- Two filter wheels
- Even better separation





Semrock quad Sedat set

# Light Sources

- Arc Lamps
  - Hg and Xe
  - Metal Halide
- LEDs
- Plasma
- Lasers
  - Generally only for collimated illumination (Confocal, TIRF)

# Arc Lamp Spectra




Hg: Brighter if your excitation spectrum matches one of the lines


Xe: More stable, longer lifetime, flat emission in visible is sometimes beneficial

# Metal Halide Arc Lamp

Exfo, Intensilight, etc.




# Metal Halide Arc Lamp





# Benefit of Liquid Light Guide

Improved Illumination Uniformity






# Plasma Lamp


LIFI; Sutter XL

- Electrodeless 'arc' source
- Very bright
- 10,000 hr life




# Plasma Lamp

LIFI; Sutter XL



# LEDs



Good in the blue / red

Not so good in the green / yellow

Long lifetime

Fast switching


# Relative lamp power

Filter Set	Excitation Filter (nm)	Dichromatic Mirror (nm)	Mercury HBO Power mW/Cm <sup>2</sup>	Xenon XBO Power mW/Cm <sup>2</sup>	Metal Halide Power mW/Cm <sup>2</sup>	LED Power mW/Cm <sup>2</sup>	Tungsten HAL Power mW/Cm <sup>2</sup>
DAPI (49) <sup>1</sup>	365/10	395 LP	23.0	5.6	14.5	0.70 (365) <sup>3</sup>	0.06 <sup>4</sup>
CFP (47) <sup>1</sup>	436/25	455 LP	79.8	25.0	76.0	26.5 (445) <sup>3</sup>	1.0
GFP/FITC (38) <sup>1</sup>	470/40	495 LP	32.8	52.8	57.5	39.2 (465) <sup>3</sup>	2.8
YFP (S-2427A) <sup>2</sup>	500/24	520 LP	20.0	35.4	26.5	10.9 (505) <sup>3</sup>	2.7
TRITC (20) <sup>1</sup>	546/12	560 LP	43.1	12.2	33.5	2.7 (535) <sup>3</sup>	1.4
TRITC (S-A-OMF) <sup>2</sup>	543/22	562 LP	76.0	31.9	67.5	6.6 (535) <sup>3</sup>	3.6
Texas Red (4040B) <sup>2</sup>	562/40	595 LP	153.7	54.4	119.5	7.9 (585) <sup>3</sup>	6.9
mCherry (64HE) <sup>1</sup>	587/25	605 LP	80.9	29.7	54.5	7.2 (585) <sup>3</sup>	4.3
Cy5 (50) <sup>1</sup>	640/30	660 LP	9.1	22.1	13.5	14.9 (635) <sup>3</sup>	4.5

<sup>1</sup>Zeiss Filters   <sup>2</sup>Semrock Filters   <sup>3</sup>LED Peak Wavelength   <sup>4</sup>Tungsten-Halogen Lamp Voltage = 12.2 V

Summary: Metal Halide best; LEDs still not quite bright enough.

# Matching Filters and Fluorophores




<http://probes.invitrogen.com/resources/spectraviewer/>


<http://fluorescence.nexus-solutions.net/frames6.htm>

<https://www.omegafilters.com/curvo2/index.php>

# Fluorescent molecules




Systems of conjugated bonds  
that share electrons



Larger system → longer wavelength

# Parameters of fluorescent molecules

- Excitation & emission maxima
- Extinction coefficient  $\epsilon$   
 *$\propto$  absorption cross section*  
 $\epsilon \approx 50,000\text{--}100,000 \text{ M}^{-1}\text{cm}^{-1}$
- Fluorescence quantum yield  $Q_f$   
*= # Photons emitted / # photons absorbed*  
 $Q_f \approx 25\text{--}90\%$   
*Brightness  $\propto \epsilon Q_f$*
- Photo-bleaching quantum yield  $Q_b$   
*= average # of photons emitted per molecule before bleaching.*  
*Depends on environment.*  
 $\propto Q_f / Q_b$




# Parameters for some common fluorophores

Dye	$\lambda_{\text{ex}}$	$\lambda_{\text{em}}$	$\varepsilon$	QY	brightness
DAPI	350	470	27000	0.58	15.7
Fluorescein	490	520	67000	0.71	47.6
Alexa 488	494	517	73000	0.6	43.8
Rhodamine	554	573	85000	0.28	23.8
Cy3	554	568	130000	0.14	18.2
Cy5	652	672	200000	0.18	36
GFP	488	507	56000	0.6	33.6
mCherry	587	610	72000	0.22	15.8
CFP	433	475	32500	0.4	13
YFP	516	529	77000	0.76	58.5



# The Enemy: *Photo-bleaching*




Decrease in emission intensity after exposure

Exciting a molecule once has a probability  $Q_b$  of killing it


Each molecule will emit only a finite number of photons

# Photo-bleaching


Photostability varies between dyes



# Photo-bleaching of fluorescent proteins



mCherry  
Single-exponential bleaching



Emerald  
Double-exponential bleaching  
Fast- and slow-bleaching populations?


# What to do about photo-bleaching?

- Select fade-resistant dyes
- Label densely
- Decrease bleaching by *anti-fade mounting media*
  - Glycerol
  - Oxygen scavengers
  - Free-radical scavengers
  - Triplet state quenchers

Note: some anti-fade agents quench some dyes.




- Budget the photons you have
  - Only expose when observing
  - Minimize exposure time & excitation power
  - Use efficient filter combinations
  - Use highly QE, low noise camera
  - Use simple light path

# Effect of mounting medium on FITC bleaching




Ono et al. 2001, *J. Histochem Cytochem.* **49**: 305-311

# Effect of mounting media on Alexa bleaching



# Fluorophore saturation



Fluorescence lifetime is  $\sim 1-5\text{ns}$

Once illumination intensity is high enough to excite the fluorophore as soon as it deexcites, further intensity increases will not increase brightness

Usually only a problem for confocal


# Factors affecting overall brightness

- Intrinsic brightness
- Spectrum of arc lamp/lasers
- Lamp/laser power
- Filter set transmission
- Quantum efficiency of detector
- Photobleaching
- Quenching / maturation / other dye-specific effects




# pH dependence of dyes


Mixed Fluorescein and Rhodamine




YFP variants




# Brief discussions of various fluors




# Traditional small molecule dyes




Alexa 350




Alexa 430




Alexa 488



Alexa 532




Alexa 546




Alexa 568

# Fluorescent dyes in Biology



FITC




Texas Red


- Protein labeling: couple to amino- or sulfhydryl groups
- Direct and indirect (immuno-) fluorescence

# Fluorescent labeling


Direct immunofluorescence:  
labeled antibodies against target




Direct labeling (& microinjection)  
of target molecules




**Indirect immunofluorescence:**  
*Unlabeled antibodies against target*  
*Labeled antibodies against those antibodies*



# DNA Probes




Ethidium Bromide  
~30 fold enhancement





DAPI  
Hoechst 33258  
Hoechst 33342  
~20 fold enhancement

# Other probes

## ER-Tracker™ Blue-White DPX



## MitoTracker Red CMXRos




Probes for Golgi, lysosomes, and peroxisomes are also available

## Small molecules – pros / cons

- 1000s available – huge spectral range
- Easy to acquire
- Precisely tailored properties, including environmental sensitivity
- Require fixing and staining, which can lead to artifacts
- Potential self-quenching and environmental sensitivity




# FIAsH/ReAsH




Example:

Newly synthesised connexins (ReAsH:Red) are added to the outer edges of existing gap junctions (FIAsH:Green). Gaietta et al 2002




# Quantum dots

- “Artificial atoms” composed of small semiconductior nanocrystals



# Quantum dots - spectra





## Quantum dots – pros / cons

- Little to no photobleaching
- Very bright
- Can use single excitation wavelength for multiple dyes
- Narrow emission spectra
- Large compared to small molecule dyes
- Single quantum dots blink
- Problems with non-specific binding

# Fluorescent Proteins and Genetically Encoded Tags

See future lecture!



# Resources

[www.microscopyu.com](http://www.microscopyu.com)

[micro.magnet.fsu.edu](http://micro.magnet.fsu.edu)

[www.chroma.com](http://www.chroma.com) (esp. their handbook on filter design)

[www.probes.com](http://www.probes.com) (esp. their handbook/catalog)



Douglas B. Murphy “Fundamentals of Light Microscopy and Electronic Imaging”

James Pawley, Ed. “Handbook of Biological Confocal Microscopy, 3rd ed.”

## Acknowledgements

Nico Stuurman / Mats Gustafsson / Mike Davidson

# Multi-photon excitation



Brad Amos, MRC, Cambridge