

The background of the slide is a high-magnification fluorescence microscopy image. It shows a dense field of cells, likely neurons or glial cells, with bright orange or yellow fluorescent structures visible within them. These structures appear to be complex, branching, or filamentous, possibly representing specific organelles or protein expressions. The overall color palette is dominated by the warm tones of the fluorescence against a dark background.

# Designing a Microscopy Experiment

Kurt Thorn, PhD  
Director, NIC@UCSF

Image from Susanne Rafelski, Marshall lab

# The Chinese-menu view of imaging

## Imaging Methods

Wide-Field

TIRF

Laser-scanning confocal

Multi-point confocal

Multi-photon confocal

## Contrasting techniques

Brightfield, phase, DIC

Immunofluorescence

Physiological dyes

Fluorescent proteins

FRAP, FLIP, photoactivation

FRET, FLIM

**Experiment:** Time Lapse, 3D, multi-point, multi-wavelength, ...

# Standard microscope capabilities

Like all rules, these were made to be broken, but only if you have fancy equipment!

- Resolution:  $\sim 200\text{nm}$  in X and Y,  $700\text{ nm}$  in Z
- Sensitivity:  $<100$  photons
- Linear detection – quantification is possible
- Video rate acquisition
- 4-5 color imaging

By far the most important part:  
*the Objective Lens*



Obviously, we care about the magnification. What other parameters are important?

# Working Distance

## Objective Working and Parfocal Distance



In general, high NA lenses have short working distances

However, extra-long working distance objectives do exist


Some examples:

10x/0.3 WD = 15.2mm


20x/0.75 WD = 1.0mm

100x/1.4 WD = 0.13mm

# Numerical Aperture



4X / 0.20 NA  
 $\alpha = 11.5^\circ$



100X / 0.95 NA  
 $\alpha = 71.8^\circ$

$$NA = n \sin(\alpha)$$

where


$\alpha$  = light gathering angle

$n$  = refractive index of sample

# Resolution of the Microscope

limited by the point-spread function


- Microscope objective collects a limited cone of light from the sample
- This limits the resolution achievable by the microscope
- Resolution can be measured by the blurring of a point object → the point-spread function



# Resolution of the Microscope


limited by the point-spread function

- Microscope objective collects a limited cone of light from the sample
- This limits the resolution achievable by the microscope
- Resolution can be measured by the blurring of a point object → the point-spread function



Y

X



Z

# Resolution of the Microscope

limited by the point-spread function

Resolution: X-Y,  $0.61\lambda / \text{NA}$ ; Z,  $\lambda n / \text{NA}^2$

Resolution for some common objectives, in nm:

<b>NA</b>	<b>X-Y</b>	<b>Z</b>
0.3	1017	16830
0.75	407	2690
0.95	321	1680
1.4	218	770

# Light-gathering power

Light-gathering power goes as the square of NA

All things being equal, a higher NA lens will give a brighter image

Increasing magnification generally decreases brightness as light is spread out over more pixels

<b>NA</b>	<b>Brightness</b>
0.3	0.09
0.75	0.56
0.95	0.90
1.4	1.96

# Choosing an objective


- Questions:
  - What resolution do you need?
  - How bright is your sample?
- For high resolution, you'll need high NA.
- For dim samples, you'll want high NA, regardless of resolution, to maximize light-gathering.
  - Dim, low-resolution samples (e.g. protein abundance in nucleus): bin camera to trade off resolution for brightness

# Choosing an objective

- Questions:
  - What resolution do you need?
  - How bright is your sample?
- When to use low NA?
  - Bright samples at low resolution / low magnification
  - If you need long working distance
  - If spherical aberration is a concern
  - If you want large depth of field to get whole structures in focus at once (avoid Z-stacks)

# NA, Z-resolution, and Z-stacks


Idea – want to record total nuclear fluorescence



# Confocal Microscopy


- Confocal microscopy has the same resolution as widefield, but eliminates out-of-focus light.
- This improves contrast for thick, heavily stained specimens.
- However, it usually comes at a cost in sensitivity.

# Fluorescence Illumination of a single point




Problem – fluorescence is emitted along entire illuminated cone, not just at focus

# The confocal microscope




What do you get?


### Axial PSF Intensity Profiles



# Confocal vs. Widefield




Confocal




Widefield

20  $\mu\text{m}$  rat intestine section recorded with 60x / 1.4NA objective

# Confocal vs. Widefield




Confocal



Widefield

Tissue culture cell with 60x / 1.4NA objective

# The confocal microscope



Scan excitation spot point-by-point to build up image

Problems:



Slow (~1 sec to acquire an image)

Low light efficiency (due to use of PMT as detector)


Solution:

Use multiple pinholes and a camera

# Spinning Disk Confocal





# Out-of-focus rejection



## Relative Sensitivity

- |                           |     |
|---------------------------|-----|
| • Widefield               | 100 |
| • Spinning-Disk Confocal  | 25  |
| • Laser-scanning Confocal | 1   |
- See Murray JM et al, J. Microscopy 2007 vol. 228 p390-405

# Multi-photon excitation



Brad Amos, MRC, Cambridge



Multi-photon excitation does not excite out-of-focus light, so you can get rid of pinhole

# Total Internal Reflection: TIRF



Thinnest optical sections:  
Images ~100nm section adjacent to coverslip

Illuminate  
through the  
objective

Sample on  
cover glass




# Epifluorescence vs. TIRF



Jaiswal et al 2002; cells loaded with FITC-dextran

# Which imaging technique should I use?




# Microscope choice

- **Epifluorescence** – routine work, low magnification, or thin samples where you don't need high-resolution 3D reconstruction
- **TIRF** – samples at the membrane or otherwise at the coverslip surface; very high signal-to-noise; single molecule imaging
- **Spinning Disk Confocal** – Live tissue culture cells, yeast, etc, or thin ( $<30\ \mu\text{m}$ ) tissue sections when you need 3D reconstructions
- **Laser-Scanning Confocal** – Thick tissues or specimens



# Sample preparation and mounting

- Mounting media serve several purposes:
  - Stabilizing the sample
  - Preventing photobleaching
  - Clearing the sample
  - Matching refractive index


# Index Mismatch & Spherical Aberration




# Index Mismatch & Spherical Aberration



$z=0\ \mu\text{m}$



$z=25\ \mu\text{m}$



$z=50\ \mu\text{m}$


# What can you do about spherical aberration?

- Use 0.17 mm coverslips ( $\sim$  #1.5)
- Work close to the coverslip
- Match lenses to the refractive index of your samples, and vice versa
  - For aqueous samples, use water immersion / water dipping lenses
  - For fixed samples and oil immersion lenses, mount your sample in a medium with  $n = 1.515$
- Adjust objective correction collar when available
- Use lower NA lenses

# Clearing

- Clearing media dissolve lipids to make samples more transparent
- Can be important for thick samples and tissues
- Commonly used:
  - BABB = 1:2 Benzyl Alcohol : Benzyl Benzoate
  - Methyl Salicylate

# Sample Preparation




Samples imaged with 20x / 0.75 air objective on spectral confocal

Sections acquired ~ 50  $\mu\text{m}$  into tissue


Embryonic mouse lungs; samples from Nan Tang, Martin Lab

# Sample Preparation


Methyl Salicylate



TDE



Vectashield



Samples imaged with 40x / 1.3 oil objective on spectral confocal

Sections acquired ~ 50  $\mu\text{m}$  into tissue

Embryonic mouse lungs; samples from Nan Tang, Martin Lab

## Dye choices – Fixed samples

- Common filter set is DAPI / FITC / Rhodamine / Cy5
- Dye choices:
  - DAPI / Hoechst / Alexa 350 / Alexa 405
  - Alexa 488
  - Rhodamine / Alexa 546 / Alexa 568
  - Cy5 / Alexa 647 / Atto 647
- More than four colors probably requires special filters or spectral imaging.

## Dye choices – Live samples


- Common filter sets: GFP / mCherry, CFP / YFP, CFP / YFP / RFP
- Two-color choice: GFP / mCherry
- Three-color: CFP / GFP / mCherry or CFP / YFP / mCherry or BFP / GFP / mCherry
- Four-color: BFP / CFP / YFP / mCherry or Sapphire / CFP / YFP / mCherry
- Five-plus colors: possible but tricky, probably requires custom filters or spectral imaging.

## Time and noise - tradeoffs


- The number of photons collected by the camera generally determines the amount of noise in your image
- Noise = square root (# of photons)
- Doubling signal to noise ratio requires 4-fold increase in exposure

# What does this look like?


With 5 e<sup>-</sup> camera read noise



1000 photons / pixel




100 photons / pixel




10 photons / pixel

# Noise and resolution

Theoretical perfect data




Two spots separated by  
diffraction limit




Slightly oversampled

# Noise and resolution


With shot noise



1000 ph/pixel at peak




100 ph/pixel at peak




10 ph/pixel at peak

# Noise and resolution


Expected error bars with shot noise



1000 ph/pixel at peak



100 ph/pixel at peak




10 ph/pixel at peak

# Noise and resolution

- High resolution and precise quantitation both require lots of light
- This means bright samples or long exposures
- This may cause problems with photobleaching and phototoxicity
- Be aware of potential tradeoffs between precision, speed, and photobleaching


# Colocalization



Utp13-GFP

Sik1-RFP

Merged



Cbf2-GFP

Sik1-RFP

Merged

Measures co-occurrence within the resolution limit of the microscope.

Does not say anything about molecular interaction

# Nothing beats good data

- Think about what data you need before you take it.
- Do you need
  - Time resolution?
  - Spatial resolution?
  - Intensity resolution?
  - Day-to-day reproducibility?
  - Spatial uniformity?
- You can fix a lot of problems with post-processing, but it's better to fix problems in the data collection!

# If you care about it, you should measure it!

- Spatial uniformity
  - Illumination and detection is not uniform over the field of view of the microscope.
  - Can be measured and corrected with a shading image.
  - Photobleaching may make this hard
- Temporal uniformity
  - Lamp power and alignment fluctuates from day to day
  - Can measure
  - But best to do experiments same day / same session

# Think about data storage

- Databases are good, but cumbersome
- Save in manufacturer's native format so metadata is preserved
- If not using a database, systematic file names and notes on sample identity are a good idea

# References

- Slides: <http://nic.ucsf.edu/edu.html>