

# Fluorescence Optics

Nico Stuurman  
Microscopy Course UCSF  
April 14, 2010

# Why fluorescence?


- High contrast
- Signal against dark background
- Highly specific, multi-color labeling
  - GFP etc.
  - Antibodies
  - Live imaging
- GFP etc.
- Quantitative
- Sensors for [Ca], pH, ...

QuickTime™ and a decompressor are needed to see this picture.

# What is it?

Sir John Frederick William Herschel, 1854: Though perfectly transparent and colorless when held between the eye and the light, or a white object, it yet exhibits in certain aspects, and under certain incidences of the light, an extremely vivid and beautiful celestial blue colour, which, from the circumstances of its occurrence, would seem to originate in those strata which the light first penetrates the liquid.....


# Excitation/Emission



Emission light is longer wavelength (lower energy)  
than excitation light


# Fluorescence

## Jablonski diagram




# Fluorescence

## Jablonski diagram




# Fluorescence Spectra



Alexa 488



# Polarization



Native GFP crystals



Shinya Inoué

# Multi-photon excitation





Brad Amos, MRC, Cambridge

# The Epifluorescence Microscope



Ploem




# Types of Filters

- Absorption (“colored”) glass
- Interference (thin-film coatings) Filters
- Acousto Optical Filters
- Liquid Crystal Filters


# Colored Glass Filters

- Cheap
- Sturdy
- Independent of angle of incidence
- Small selection
- Spectra have poor slope and poor peak performance
- Autofluorescence
- Absorb → Get Hot


# Interference Filters



# Interference Filter Design (multiple cavities)




# Highly sensitive to angle




Positive spin on this pitfall:  
Tunable Bandpass Filters!

(Semrock *VersaChrome*)


Bandpass




Longpass




%T



OD



# Matching Filters and Fluorophores




<http://probes.invitrogen.com/resources/spectraviewer/>



<http://fluorescence.nexus-solutions.net/frames6.htm>

<https://www.omegafilters.com/curvo2/index.php>


## Choose filters that separate fluorophores




Two different UV filter sets



Choose filters that maximize excitation and emission




Choose filters that maximize excitation and emission



# Newer hard-coatings are great!



# are great!



# Blocking


Interference filters have finite stop bands

Unblocked bandpass interference filter




To block unwanted transmission from UV to IR, filter makers add absorption glass to the filter.


Often excitation filters are blocked,  
but emission filters unblocked.  
→ Red autofluorescence or room light  
may get through your blue emission filter



# Faster Wavelength Selection: Multiple Band Pass Filters & Filter



# Faster Wavelength Selection: Multiple Band Pass Filters & Filter



# Filter schemes

## Single wavelength sets

- Most efficient
- Best separation
- Very slow to change  $\lambda$



Transmission



## Multi-band filters

- Multi-band everything
- See all colors at once
- For color cameras
- Bad crosstalk
- “Pinkel” scheme

Multi-band dichroic

Multi-band emitter

Single- $\lambda$  excitors

- Excitation filter wheel
- Separate image at each wavelength
- Better separation

• “Sedat” scheme

Multi-band dichroic


single-band emitters

Single- $\lambda$  excitors



- Two filter wheels
- Even better separation



Wavelength




Chroma triple Pinkel set



Semrock quad Sedat set

# Acoustical Optical Tunable Filter




- Switches and modulates intensity
- Fast! (sub-microseconds)
- Mainly used for excitation laser light
- Polarization dependent


Also: AOM or Bragg cell

# Optical Filter

Figure 2





# Liquid Crystal Filters




- Example: Lyot filter:  
Uses Birefringence  
and polarizers
- Shifts in (100) ms time  
range
- Maximum transmission  
is 50%, blocking max  
 $10^{-5}$

# Light source spectra

Mercury Arc Lamp UV and Visible Emission Spectrum




Laser Illumination Source Emission Spectra



Solid-state lasers: many, many lines available

LEDs are here!

# Koehler illumination



# Thanks!

- Mats Gustafsson
- Kurt Thorn
- Jennifer Waters
- <http://micro.magnet.fsu.edu/>
- <http://www.microscopyu.com>
- <http://olympusmicro.com>
- <http://zeiss-campus.magnet.fsu.edu/>
- <http://www.chroma.com> (Filter Handbook!)

