

C-Trap data analysis

An abstract, stylized graphic on the right side of the slide. It features a complex, branching, and somewhat circular structure in shades of red, pink, and white. The structure has a textured, almost pixelated or splattered appearance. It is set against a solid blue background that covers the entire slide.

An introduction to data analysis using Lakeview
and Pylake

Aafke van den Berg, PhD
Sr Customer Success Data Scientist




Table of Content


- Recourses for data analysis:
 - Lakeview
 - Pylake
 - Harbor
- Pylake example script
- Lakeview demo

Data analysis


Analyzing C-Trap hdf5 data using **Lake**view & Py**lake**

Exporting data from Bluelake

Data Visualization and management




Instrument Control



- Each time Bluelake is opened, a new **session** is created that contains **all data streams**.
- Note: due to the size of session files, they cannot be kept permanently.
- Sessions can only be opened in **Session Viewer**, which can be found on the desktop of the Controller PC.

Selected channels in marker can be exported as .h5 file

Marked data




- Visualize/show




- Download/include channel when exporting file



List of data that can be exported



- List of Scans, Kymographs, Point scans, FD curves, etc that were collected during the experiment
- The data can be exported as .h5 file.
- Scans and kymographs can also be saved as an image

hdf5 files




- Data exported from Bluelake is stored in the hdf5 (.h5) format .
- The HDF5 files contain groups (Force HF, Distance,...) that can hold datasets (Force 1x, Distance 1, ...)
- More information about HDF5 files can be found at <https://support.hdfgroup.org/HDF5/doc/H5.intro.html>
- HDF5 files can be opened with various data analysis software, such as Matlab, Python, origin, prism, ...
- LUMICKS has developed
 -  **Lakeview** for direct visualization of hdf5 files
 -  **PyLake** for data analysis

C-Trap data analysis

Lakeview


- Quickly visualize C-Trap h5 data
- Export data to a different format
- Perform Kymograph analysis



- Download Lakeview [here](#)

PyLake


- Python package for visualizing and analysing C-Trap h5 data, for example:
 - Analyzing fd curves
 - Computing the Piezo distance
 - Example analysis workflows



- [PyLake documentation](#)

harbor

- Script-sharing platform for C-Trap users
- Download data analysis and automation scripts




- harbor.lumicks.com

Example data output from Pylake


Data channel over time


Example: sudden distance changes over time reflects the unwrapping of nucleosomes





Data channel vs. another data channel

Example: force-distance relationship of dsDNA, analyze DNA compaction








Force-extension curves



Stretching DNA




- When performing force-extension measurements, you often want to know
 - What is the length of my double stranded DNA?
 - What are the mechanical properties of my DNA?
 - Level of compaction
 - ...
- This information can be obtained by fitting a model to the force extension curve




The worm-like chain model


- The Worm-like chain (WLC) model is a good model for fitting DNA force-extension curves
- Variants of the WLC model:
 - WLC (up to 5-10 pN)
 - Extensible WLC (up to 30 pN) -> **most common**
 - Twistable WLC (up to 50 pN)
 - ...
- Choose model with as few parameters as possible



Extensible WLC



Extensible worm-like chain (up to ~ 30 pN)
Odijk 1995



In Pylake:


`ewlc_odijk_force`

`ewlc_odijk_distance`


$$d(f) = Lc_{DNA} \left(1 - \frac{1}{2} \sqrt{\frac{kT}{fLp_{DNA}}} + \frac{f}{St_{DNA}} \right)$$

Parameters of the extensible WLC: Contour length

- The contour length, L_c , is the length of the DNA




- The contour length can be estimated from the FD curve




Parameters of the extensible WLC: Persistence length, L_p

- Measure of stiffness: 'the persistence length is the typical length along which the chain forgets its previous orientation.'
- Typically 50 nm for double-stranded DNA
- Value depends on salt concentration (Baumann 1997)




Parameters of the extensible WLC: Stretch modulus

- Higher St means more resistance against stretching
- St important for higher forces
- Typical value for dsDNA: 1500 pN



Extra parameter: Force offset

- The force offset can be non-zero when the force was not reset before recording the FD curve
- The force offset vertically shifts the fd curve
- Note: The force offset can be negative, this can give weird effects when using the absolute force. Therefore, it is better to use Force 2x






Nucleosome unwrapping

DNA tethering kit



Nucleosome array



Concluding remarks

- Pylake can be cited using the following DOI: 10.5281/zenodo.4280789
- Both questions and feedback can go to the support website or support@lumicks.com

A large, diverse group of people, likely a company team, posing on a set of stairs. They are all smiling and making various hand gestures, including peace signs, thumbs up, and open palms. The group is multi-ethnic and multi-gender. The background shows a modern building with large windows and a staircase railing. The overall mood is energetic and positive.

We are LUMICKS!

A team of experts, passionate to help you