



Microscopy Optical Sectioning



Stephen Ross Ph.D.
UCSF Principles and Practice
Of Light Microscopy
May 3rd, 2010

The Eyes of Science



- **Confocal Microscopy In Principle**
- **Point Scanning Confocal**
- **Spinning Disk Confocal**
- **Swept Field Confocal**
- **Total Internal Reflection Fluorescence (TIRF)**

Broad definitions –An introduction to understanding the terminology

1. Define optical section in terms of the mathematics:

- Lateral resolution: for diffraction limited techniques: where d is the smallest resolvable unit, λ is the wavelength of light, and NA is the numerical Aperture of the objective lens.

$$d_{xy} = \frac{1.22\lambda}{(2NA)}$$

Example: For a 1.4NA objective at 550nm, the resolution limit is 240nm.


Broad definitions –An introduction to understanding the terminology

1. Define optical section in terms of the mathematics:

- Axial resolution (focal volume): Again, for diffraction limited techniques, the axial resolution is given as: where η is the refractive index of the media.


$$d_z = \frac{2\lambda \cdot \eta}{(NA_{obj})^2}$$

Example: For a 1.4NA objective in oil/oil, at 550nm, the resolution limit is 850nm.





Point Spread Function

- PSF defines the propagation of electromagnetic radiation or other imaging light rays from a point source or point object.
- The degree of spreading of the point object is a measure of the imaging system.
- The effect of the imaging system on a point object is termed the convolution of the optical system.
- Displayed by imaging a sub micron point in 3D (Z stack) and displaying as a volumetric airy disc. (X-Z projection)




Diffraction limited
Point Source



Optical Convolution of
Point Source

Multiple Focal Planes




The Eyes of Science



Conventional light microscopy is limited by multiple focal planes in the same section, therefore:

- Physically cut thinner sections
 - Technically difficult
 - Highest resolution
 - Highest sensitivity
- Optical sectioning using confocal microscopy
 - Technically simple
 - Multiple section planes available
 - Complex computing

Principal of Confocal microscopy



The Eyes of Science




Pinhole size and resolution

- Pinhole size should be related to width of point spread function
- Width of point spread function = resolution of lens \times magnification of lens = 1 Airy unit
 - 100x / 1.4 NA: resolution = 220nm, so 1 Airy unit = 22 μm
 - 40x / 1.3 NA: resolution = 235nm, so 1 Airy unit = 9.4 μm
 - 20x / 0.75 NA: resolution = 407nm, so 1 Airy unit = 8.1 μm
 - 10x / 0.45 NA: resolution = 678nm, so 1 Airy unit = 6.8 μm


Confocal and Widefield Fluorescence Microscopy




(a)




(b)




(c)



(d)



(e)



(f)

Figure 1

Diffraction limited volume :


$$NA = 1.4$$

$$\eta \text{ (Refractive Index)} = 1.515$$

$$\lambda \text{ (Wavelength)} = 550 \text{ nm}$$

$$d_{xy} = \frac{1.22\lambda}{(2NA)} = \frac{1.22 \cdot 0.55 \mu m}{(2 \cdot 1.4)} = 0.24 \mu m = 240 \text{ nm}$$

$$d_z = \frac{2\lambda \cdot \eta}{(NA_{obj})^2}$$
$$= \frac{2 \cdot (0.55 \mu m) \cdot 1.515}{(1.4)^2} = 0.85 \mu m = 850 \text{ nm}$$



Components of a confocal




Light source
Generally laser based



Detector
PMT or CCD




Controller and
image acquisition
computer




Scan Head
w/pinholes


Microscope

3D reconstruction



Once the “Z-stack” has been collected and sent to the computer, Sections of the stack can be reconstructed into a 3D volume.






The Eyes of Science




- **Confocal Microscopy In Principle**
- **Point Scanning Confocal**
- **Spinning Disk Confocal**
- **Swept Field Confocal**
- **Total Internal Reflection Fluorescence (TIRF)**

How a Point Scanner Works



The pinhole will only provide information about a single point in space. Therefore move specimen or move the exciting light




X and Y Galvo Mirrors scan the whole field in a raster pattern

- Note that the mirrors also descan the returning light
- Many systems have been developed to control the scan mirrors
- Also motorize Z-axis to sample points in vertical axis


Commonly used Lasers

- Argon Ar 457, 488, 514 nm
- Krypton-Ar Kr-Ar 488, 568, 647 nm
- Helium-Neon He-Ne lots single line sources
543, 594, 633 nm
- He-Cadmium He-Cd 325 - 441 nm
- Diodes 350, 405, 442, 473, 488, 491, 514,
532, 561, 593, 638, more coming!

Photomultipliers (PMTs)



- Maximize the light budget returning from the scan head by intensification.
- Subject to electrical noise to much gain will degrade image and demand excessive averaging






PMTs must be matched with detection wavelength

The Eyes of Science



Computer builds digital image from scanned data




The Eyes of Science



Due to PMT Noise, Averaging is commonly used to enhance S/N





1 frame



8 frames

Multi-photon excitation



Brad Amos, MRC, Cambridge



Multi-photon excitation does not excite out-of-focus light, so you can get rid of pinhole

The Eyes of Science



- **Confocal Microscopy In Principle**
- **Point Scanning Confocal**
- **Spinning Disk Confocal**
- **Swept Field Confocal**
- **Total Internal Reflection Fluorescence (TIRF)**

Spinning Disc Confocal



- **Confocal Microscopy In Principle**
- **Point Scanning Confocal**
- **Spinning Disk Confocal**
- **Swept Field Confocal**
- **Total Internal Reflection Fluorescence (TIRF)**

SFC Design Overview


- Supravideo rate, high frequency **slit** scanner.
- Variable pinhole, linear array **field** scanner.
- CCD camera-based detector.
- Linear scanning at high frequency minimizes phototoxicity.



SFC Design Overview, cont'd.

- Variable pinhole sizes match various objective lenses.
- Linear pinhole array (1-D) and separate excitation and emission light path minimize emission crosstalk.

SFC Aperture Positions






Emission slits
and pinholes

Excitation slits
and pinholes


The Eyes of Science




Laser beam is shaped into a line, then scanned




Slit scanner lightpath (simplified)






The Eyes of Science







Low frequency
galvo sweep




High frequency
galvo sweep

Field Scanner lightpath (simplified)






Field-scanner lightpath




The Eyes of Science





Galvo "Sweep"





The Eyes of Science




Piezo movement of Pinhole image

CCD Frame #1 (or odd #s)




CCD Frame #2 (or even #s)





Piezos move a fixed distance

Time to move depends on exposure time

It takes 2 frames for one complete "piezo sweep"





Galvo sweeps **multiple** cycles/frame,
and piezos sweep only **once** per frame



Galvo cycles multiple times / frame

The Eyes of Science






Field-scanning combines multiple galvo sweeps
with bidirectional piezo movement

The Eyes of Science




← Piezo step →




Galvo "Sweep" ↑
↓

The Eyes of Science







Galvo, Piezos stationary



Sweeping galvo only



Sweeping galvo & piezos




Piezo sweep sync
with frame rate

The Eyes of Science




Performance


- Slit scan imaging to 1000Hz with appropriate detector.
- Pinhole scan imaging to 100Hz.
- Axial resolution from $\sim 0.4\mu\text{m}$ to $1.2\mu\text{m}$ over 6 scanning modes.



10 minute
timelapse (35um
slit) Beta cell with
eGFP labeled
insulin.




Calcium sparks at
80 FPS in smooth
muscle cells
(Fluo4)



The Eyes of Science





Drosophila Embryo
Heart Beating 50
FPS, GFP
expressed in
smooth muscle.

The Eyes of Science




- **Confocal Microscopy In Principle**
- **Point Scanning Confocal**
- **Spinning Disk Confocal**
- **Swept Field Confocal**
- **Total Internal Reflection Fluorescence (TIRF)**

What is Evanescent Wave or Total Internal Reflection Fluorescence (TIRF) Microscopy?

- An optical phenomenon that can occur when light strikes the interface between two media of different refractive indices
- When light is totally internally reflected it does not propagate through the sample, but causes an electromagnetic field (termed the evanescent wave) that extends in the z-direction
- The evanescent wave intensity decreases exponentially, and extends only a few hundred nanometers into the second media


Glass ($n=1.515$)



Live cells ($n=1.33-1.38$)

Criterion for TIR

- $n_1 > n_2$
- The angle of incidence must be \geq the critical angle (θ_c)



$$E_z = E_0 \exp(-z/d_p)$$

$$d_p = \frac{\lambda_i}{2\pi n \sqrt{\sin^2 \Theta_i - (n_2/n_1)^2}}$$

- The energy of the evanescent wave at any point in z is a function of location in the field and the penetration depth
- The penetration depth of the evanescent wave is a function of the illumination wavelength.
- only 50nm to 100nm is typically useful illumination

Why is this useful for fluorescence microscopy?


- The evanescent wave will typically only excite fluorophores within 50nm-80nm of the cover slip.
- Excitation of this small proportion of fluorophores in the thin optical section results in an extremely large increase in the S/N
- As a result of the increase in S/N very small signals, such as single molecule fluorescence, can be resolved above the noise floor

TIRF is all about

SIGNAL to NOISE!

The Eyes of Science






SIGNAL to NOISE
Epi-fluorescence
=1.3
For this specimen


Epi-Fl Image

5mm deep solution of 200nm fluorophor coated polystyrene micro-spheres.

3D Intensity Histogram of image on the left



Low Intensity (black)  High Intensity (white)





TIRF Image

5mm deep solution of 200nm fluorophor coated polystyrene micro-spheres.

**SIGNAL to NOISE
TIRF
=35
For the same specimen**

3D Intensity Histogram of image on the left




Low Intensity (black)  High Intensity (yellow)

What are some of the applications of TIRF ?

- Endocytosis and Exocytosis
- Dynamics of membrane associated proteins
- Protein arrangement at membranes
- Focal adhesions
- Growth cone migration
- Receptor/Ligand interactions
- Biophysical studies of single molecule behavior
- ???? For you to figure out

Objective lens TIRF Method




High NA objective lenses are necessary for Total Internal Reflection, utilizing the lens method, so that light can come out of the lens at angles at or beyond the critical angle for TIR

$$\text{NA} = n \sin \theta$$

Numerical
Aperture

Mounting
Media
Refractive
Index


Angle of the
Cone of
illumination




The Eyes of Science



High Numerical Aperture Objective TIRF





**Special objective lenses provide the NA
Necessary to achieve Total Internal Reflection
Through the objective lens.**



The Eyes of Science






The Eyes of Science





Ti-SFL



COS Cell Expressing GFP Clathrin

(Taken with "white light" TIRF)



Images by: Prof. Dan Axelrod, Univ. of Michigan


The Eyes of Science



Aplysia Growth Cone, primary culture

Alexa 488 Actin

Camera: Roper Coolsnap HQ, 0.1Hz



Imaged by Andy Schaefer, Paul Forscher Lab. , Yale University

The Eyes of Science




Aplysia Growth Cone, primary culture

Alexa 488 Actin

Alexa 568 Tubulin

Camera: Roper Coolsnap HQ, 0.1Hz




Imaged by Andy Schaefer, Paul Forscher Lab. , Yale University

The Eyes of Science



Aplysia Growth Cone, primary culture
Alexa 488 Actin (speckled concentration)
Alexa 568 Tubulin
Camera: Roper Coolsnap HQ, 0.1Hz



Imaged by Andy Schaefer, Paul Forscher Lab. , Yale University

The Eyes of Science



Thank You

The Eyes of Science

