


A fluorescence microscopy image of plant tissue, likely a leaf cross-section. The image shows a network of green, thread-like structures (cell walls or chloroplasts) forming a honeycomb-like pattern. Within these green structures are numerous small, bright red spots, which could be chloroplasts or other cellular components. The overall image has a dark background, making the green and red structures stand out.

# Optical Sectioning and Confocal Microscopy


Kurt Thorn

The goal: build 3D images of biological samples




Embryonic mouse lung; Nan Tang, Martin lab


# Optical Sectioning and 3D reconstruction




# Optical Sectioning and 3D reconstruction



# Optical Sectioning and 3D reconstruction




# A series of optical sections



Yeast mitochondria; Susanne Rafelski

# Fluorescence Illumination of a single point



Problem – fluorescence is emitted along entire illuminated cone, not just at focus

# Widefield fluorescence imaging






# The Point Spread Function (PSF)


The image of a point object

2D PSF  
for different defocus




3D PSF

Calculated



Measured




# The Problem of Out-of-Focus Light

## Three Options


- Prevent it from reaching the detector
  - Confocal microscopy
- Prevent it from being excited
  - Two-photon, light sheet microscopy, TIRF
- Remove it after the fact (computationally)
  - Deconvolution

# Fluorescence Illumination of a single point



Problem – fluorescence is emitted along entire illuminated cone, not just at focus


# The confocal microscope



Scan excitation spot point-by-point to build up image


# How big should your pinhole be?

Resolution is limited by the point-spread function




Y

X



Z




Airy disk radius

$$\approx 0.61 \lambda / NA$$


# How big should your pinhole be?

Want pinhole to pass entire Airy disk



Y

X



Airy disk diameter  
 $\approx 1.22 \lambda / NA$

Width of point spread function at pinhole:  
Airy disk diameter  $\times$  magnification of lens

## How big should your pinhole be?


- Width of point spread function at pinhole =  
Airy disk diameter  $\times$  magnification of lens = 1 Airy unit  
= resolution of lens  $\times$  magnification of lens  $\times$  2
  - 100x / 1.4 NA: resolution = 220nm, so 1 Airy unit = 44  $\mu\text{m}$
  - 40x / 1.3 NA: resolution = 235nm, so 1 Airy unit = 19  $\mu\text{m}$
  - 20x / 0.75 NA: resolution = 407nm, so 1 Airy unit = 16  $\mu\text{m}$
  - 10x / 0.45 NA: resolution = 678nm, so 1 Airy unit = 14  $\mu\text{m}$

# Light sources

Excitation light must be focused to a diffraction limited spot


Could be done with an arc lamp and pinhole – but very inefficient

Enter the laser:  
Perfectly collimated and  
high power






# Scanning



Changing entrance angle of illumination moves illumination spot on sample


The emission spot moves, so we have to make sure pinhole is coincident with it

# Confocal optical path




What do you get?

### Axial PSF Intensity Profiles



# Detectors - PMTs

- Must be fast – confocal beam spends only a few  $\mu\text{s}$  on each pixel
  - Photomultiplier tubes



Pulse width for  
single photon  
 $\sim 10\text{-}100\text{ns}$

Very linear


Very high gain  
 $\sim 0$  read noise

# Detectors - PMTs


- Peak QE for standard PMT: ~30%
- Front-illuminated CCD: ~60%
- Back-illuminated CCD: ~95%
- GaAsP PMTs: ~40-50%

However, GaAsP PMTs can be damaged by exposure to too much light.


# Detectors - APDs



- Photons produce photoelectrons
- Photoelectrons are amplified by collisional ejection of electrons
- Small dynamic range




# Hybrid Detectors




- GaAsP photocathode bombards an APD
- Very high gain
- Excellent for single photon counting
- Damaged by excessive light

# Confocal vs. Widefield



Confocal




Widefield


Tissue culture cell with 60x / 1.4NA objective



# Confocal vs. Widefield



Confocal




Widefield

20  $\mu\text{m}$  rat intestine section recorded with 60x / 1.4NA objective

# Doing more with your photons


## Spectral Detection



Prism + Adjustable Slits: Leica SP microscopes

# Doing more with your photons


## Spectral Detection




Prism + PMT-array: Zeiss and Nikon

# Doing more with your photons

## Time-Correlated Single Photon Counting



# The confocal microscope: Drawbacks



Scan excitation spot point-by-point to build up image

Problems:

Slow (~1 sec to acquire an image)

Low light efficiency (due to use of PMT as detector)

Solution:

Use multiple pinholes and a camera

# A Solution: Spinning Disk Confocal




Figure 5






Image with many pinholes at once, so fast

Use CCD as detector, so much higher QE

# The downside to the spinning disk




# The downside to the spinning disk






# The downside to the spinning disk



# The downside to the spinning disk


Limited out of focus rejection




# Pros/Cons of spinning disk

- Fast – multiple points are illuminated at once
- Photon efficient – high QE of CCD
- Gentler on live samples – usually lower laser power
- Fixed pinhole size (usually)
- Small field of view (usually)
- Crosstalk through adjacent pinholes limits sample thickness

# Out-of-focus rejection



# Examples



*Drosophila* S2 cell expressing  
GFP-H2B and mCherry-tubulin  
(Nico Stuurman and Ron Vale)



*S. cerevisiae* expressing a  
mitochondrially targeted RFP,  
Susanne Rafelski, Marshall lab

# When to use confocal?

- Confocal is not a magic bullet
  - It is extremely wasteful of photons
  - Laser-scanning confocal is 100 – 200-fold less sensitive than widefield
  - Spinning-disk confocal is ~4-fold less sensitive than widefield

# When to use confocal?

- Confocal is not a magic bullet
  - It is extremely wasteful of photons
  - High laser power generally result in more photobleaching and photodamage.
- For thin specimens, widefield epifluorescence is better – especially with deconvolution
- Confocal excels with thick, heavily stained specimens

# When to use confocal?


- What is thick?
  - A good rule of thumb is  $10 \times$  the depth of field of the objective
  - 100x / 1.4 NA: d.o.f.  $0.66\mu\text{m}$
  - 20x / 0.75 NA: d.o.f.  $2.3\mu\text{m}$
- Sample preparation is KEY for imaging thick specimens
  - Confocal does not fix scattering, refractive index mismatch, or everything else that can go wrong – it only removes out of focus light



# Sample preparation

- For fixed samples: match refractive index of mounting media to immersion oil.
  - Mount in immersion oil itself, BABB, benzyl alcohol/glycerol, 2,2'-thiodiethanol, or other high-RI mounting medium
- Clearing to remove lipids and other scattering substances is also important
- For live samples, use water immersion lenses

# Sample Preparation




Samples imaged with 20x / 0.75 air objective on spectral confocal

Sections acquired ~ 50  $\mu\text{m}$  into tissue


Embryonic mouse lungs; samples from Nan Tang, Martin Lab

# Sample Preparation


Methyl Salicylate



TDE



Vectashield



Samples imaged with 40x / 1.3 oil objective on spectral confocal

Sections acquired ~ 50  $\mu\text{m}$  into tissue



Embryonic mouse lungs; samples from Nan Tang, Martin Lab

# Clearing and mounting summary

- Both clearing and refractive index matching are important.
- BABB and methyl salicylate clear very well and give best image depth, but may disrupt cell morphology
- TDE preserves cell morphology and also allows for decent imaging depth

# Newer Clearing Techniques



## Scale



ScaleA2: 4M urea, 10% glycerol, 0.1% Triton X-100

# Newer Clearing Techniques

## CLARITY



# Resources

<http://www.microscopyu.com>

<http://micro.magnet.fsu.edu>

James Pawley, Ed. “Handbook of Biological Confocal Microscopy, 3rd ed.”

[http://www.hamamatsu.com/resources/pdf/etd/PMT\\_handbook\\_v3aE.pdf](http://www.hamamatsu.com/resources/pdf/etd/PMT_handbook_v3aE.pdf)

# Acknowledgements

Steve Ross, Mats Gustafsson