

Principles & Practice of Light Microscopy 6

Special Techniques

(FRET, FRAP , FLIP , FLIM , FCS, molecular sensors...)



(Image: T. Wittman, Scripps)

Monitoring molecular movement by microscopy

- General idea: optically mark one region of the cell, follow where it goes at later times.
- Implementation
 - Uncaging (old)
 - FRAP, FLIP: photobleaching
 - Photoactivation


FRAP

Fluorescence Recovery after Photo-bleaching



Need: to probe transport

Idea: bleach in one area,
watch recovery by transport from other areas

Fluorescence Recovery After Photobleaching (FRAP) with Green Fluorescent Protein



Measuring Cdc42 diffusion constant in yeast



Result: $d_f = (0.036 \pm 0.017) \mu\text{m}^2/\text{s}$


Marco et al. 2007 Cell 129:411-422

FLIP

Fluorescence Loss in Photo-bleaching

Need: probe connectivity

Idea: bleach in one compartment,
watch loss in connected compartments by exchange



Bleach one area repeatedly. Entire ER dims.

⇒ ER is contiguous


Photo-Activation

(Better?) FRAP/FLIP alternative

Some fluorophores can be activated by light

- Photo-uncagable dyes
- GFP-family proteins

Look for weak light
against dark background
Instead of slight dimming
of bright background



Activate a small area
Watch fluorescence spread

Photo-Activatable Proteins

Off-On


- PA-GFP, PS-CFP

Color change


- Kaede, KikGR, Eos,
- Dendra (activatable by blue)

Reversibly Switchable

- asCP, KFP (tetrameric)
- Dronpa






Dronpa – photoswitchable on and off



D R O N P A


Tracking actin flow with Dronpa



Hardware for photoactivation / photobleaching



- Need: controlled illumination of a region in the sample

Laser-scanning confocal microscope



Use mirrors to adjust angle at back focal plane – changes location of laser spot at sample


Laser-scanning confocal microscope




Laser-scanning confocal microscope

- Pros
 - can bleach / activate any shape ROI by scanning point by point
- Cons
 - Slow
 - Need a laser scanning confocal

Can control field of view with field iris



Better idea: Spatial Light Modulator




Can project any image onto sample


Digital Micromirror Device



Array of rotatable mirrors



Spatial Light Modulator



Pros: Fast, can project any arbitrary image


Cons: Intensity is relatively low

FLIM

Fluorescence Lifetime Imaging

Measure the lifetime of the excited state

- Separate fluorophores with similar spectra
- Detect environmental parameters that affect lifetime (FRET, pH, ...)



FLIM



Measurement approaches

- Frequency domain

- Modulated excitation
- Lock-in detect emission phase


- Time domain (pulsed exc.)

- Gated intensifier
 - Photon inefficient
- Time-correlated single photon counting
 - Very efficient
 - \leq one photon per pulse \Rightarrow slow



FLIM



Examples





Hepatocyte membrane-stained
with NBD, which has a
hydrophobicity-dependent lifetime
(TCSPC, 3 minutes for 300x300 pixels)

Environment-sensitive fluorophores



Many: pH, ions (Ca^{2+} , Mg^{2+} , Na^+ , K^+ , etc.), voltage, hydrophobicity,...



Fura-2








Fluorescein




MQAE

Environment-sensitive fluorophores





Imaging Cdc42 activation





Fuse sensor to GFP; I-SO/GFP ratio is proportional to degree bound

Monitoring Cdc42 activation



Fluorescence Resonance Energy Transfer




Sensitive to distances on the order of 10 nm,
20-fold below resolution limit

Distance dependence of FRET

$$E = \frac{1}{1 + (r^6/R_0^6)}$$

$$R_0^6 \propto k^2 n^{-4} Q_D J(I)$$

Overlap between donor
 emission and acceptor
 excitation
 Donor quantum yield
 Refractive index
 Orientation between
 fluorophores




For CFP-YFP,
50% transfer at $R_0 = 4.9$ nm

Good FRET pairs

- CFP/YFP – use A206R mutants if dimerization is problematic
- GFP/mCherry or other FP pairs – not so well validated
- Fluorescein/Rhodamine
- Cy3/Cy5 or Rhodamine/Cy5
- Many other small molecule pairs


FRET Theory

- $k^2 = (\cos q_T - 3 \cos q_D \cos q_A)^2$
- For rapidly tumbling molecules, can average over all possible orientations to give $k^2 = 2/3$
- But rotational correlation time for GFP is ~ 16 ns; fluorescence lifetime is ~ 3 ns




Effects of FRET

- Donor lifetime shortened
- Acceptor emission depolarized
- Donor fluorescence quenched
- Acceptor fluorescence enhanced on donor excitation




Measuring FRET

- Donor lifetime shortened
- Can measure by fluorescence lifetime imaging, but requires specialized instrumentation




Measuring FRET

- Acceptor emission depolarized
- Can measure by fluorescence polarization microscopy




Measuring FRET

- Donor fluorescence quenched
- Acceptor fluorescence enhanced on donor excitation
- Can measure by donor recovery after acceptor photobleaching
 - Easy, but very sensitive to degree of photobleaching




Measuring FRET




- Donor fluorescence quenched
- Acceptor fluorescence enhanced on donor excitation
- Can measure by quantitative measurement of acceptor enhancement on donor excitation
 - Easy, but crosstalk needs to be corrected for

A problem: crosstalk into FRET channel




Correct using measurements from CFP- and YFP- only cells


Crosstalk correction



Types of FRET experiments




Intramolecular



Intermolecular

Types of FRET experiments




Intramolecular

For intramolecular FRET, CFP and YFP are always present in a 1:1 ratio

Ratiometric imaging can be used as a rough measure of the amount of energy transfer

Types of FRET experiments





Intermolecular

For intermolecular FRET, the relative abundance of CFP and YFP is not controlled and can change over time.

Ratiometric imaging is no longer possible, and additional corrections are necessary.

Using FRET to monitor Rac activation



FCS

Fluorescence Correlation Spectroscopy


Small volume \rightarrow only a few molecules \rightarrow random fluctuations

- Study the *noise*
- Conclude about *random processes* at different *time scales*

Small excitation volume

Random molecular processes:

- Diffuse in and out
- Adopt different states
- Bind or react
- Photobleach
- ...



Study the auto-
(or cross-) correlation




Image Correlation Spectroscopy

Image FCS

Sense the random fluctuations of fluorescence within an *image* or *image sequence*

Much slower than point FCS, but get whole area

Can see where you are & deal with motion

Can do spatio-temporal analysis

Image series

Measure time and space variations in intensity







Image Correlation Spectroscopy




Areas
analyzed

Flow
vector

$\varnothing \approx 10\mu\text{m}$
diffusion
length



 Stationary fraction
 Flowing fraction
 Diffusing fraction





Spatio-temporal
autocorrelation


Optical Tweezers

Mechanically manipulate
the specimen with light

Why?



Measuring force and displacement
of a single polymerase molecule



DNA recoil after stretching

Optical Tweezers


How does it work?

Photons carry momentum

Changing photon direction requires a force

More light refracted away from high intensity region

→ Force toward intensity maximum



Alternative way to look at it:



Field energy is less in higher index

System energy depends on particle position

$$dE/dx = \text{Force}$$

Holographic Optical Tweezers

- Many traps at once
- Independently movable
- Made using a computer-controlled spatial light modulator in a pupil plane



Further reading

www.microscopyu.com

micro.magnet.fsu.edu

Molecular Probes Handbook (probes.com)

James Pawley, Ed. "Handbook of Biological Confocal Microscopy, 3rd ed."

Acknowledgements

Mats Gustafsson