


A fluorescence microscopy image of plant tissue, likely a leaf cross-section. The image shows a network of green, thread-like structures (cell walls or chloroplasts) forming a honeycomb-like pattern. Within these green structures are numerous small, bright red clusters, which could be stomata or specific cellular components. The overall image has a dark background, making the green and red structures stand out.

Optical Sectioning and Confocal Microscopy


Kurt Thorn

The goal: build 3D images of biological samples




Embryonic mouse lung; Nan Tang, Martin lab


Optical Sectioning and 3D reconstruction




Optical Sectioning and 3D reconstruction



Optical Sectioning and 3D reconstruction




A series of optical sections



Yeast mitochondria; Susanne Rafelski

Fluorescence Illumination of a single point



Problem – fluorescence is emitted along entire illuminated cone, not just at focus


Widefield fluorescence imaging



The Point Spread Function (PSF)


The image of a point object

2D PSF
for different defocus




3D PSF

Calculated



Measured




The Problem of Out-of-Focus Light

Three Options


- Prevent it from reaching the detector
 - Confocal microscopy
- Prevent it from being excited
 - Two-photon, light sheet microscopy, TIRF
- Remove it after the fact (computationally)
 - Deconvolution

Fluorescence Illumination of a single point



Problem – fluorescence is emitted along entire illuminated cone, not just at focus


The confocal microscope



Scan excitation spot point-by-point to build up image


How big should your pinhole be?

Resolution is limited by the point-spread function




Y

X



Z




Airy disk radius

$$\approx 0.61 \lambda / NA$$


How big should your pinhole be?

Want pinhole to pass entire Airy disk



Y

X



Airy disk diameter
 $\approx 1.22 \lambda / NA$

Width of point spread function at pinhole:
Airy disk diameter \times magnification of lens

How big should your pinhole be?


- Width of point spread function at pinhole =
Airy disk diameter \times magnification of lens = 1 Airy unit
= resolution of lens \times magnification of lens \times 2
 - 100x / 1.4 NA: resolution = 220nm, so 1 Airy unit = 44 μm
 - 40x / 1.3 NA: resolution = 235nm, so 1 Airy unit = 19 μm
 - 20x / 0.75 NA: resolution = 407nm, so 1 Airy unit = 16 μm
 - 10x / 0.45 NA: resolution = 678nm, so 1 Airy unit = 14 μm

Light sources


Excitation light must be focused to a diffraction limited spot

Could be done with an arc lamp and pinhole – but very inefficient

Enter the laser:
Perfectly collimated and
high power




Scanning



Changing entrance angle of illumination moves illumination spot on sample


The emission spot moves, so we have to make sure pinhole is coincident with it

Confocal optical path




What do you get?

Axial PSF Intensity Profiles



Detectors - PMTs

- Must be fast – confocal beam spends only a few μs on each pixel
 - Photomultiplier tubes



Pulse width for
single photon
 $\sim 10\text{-}100\text{ns}$

Very linear

Very high gain


~ 0 read noise

Detectors - PMTs


- Peak QE for standard PMT: ~30%
- Front-illuminated CCD: ~60%
- Back-illuminated CCD: ~95%
- GaAsP PMTs: ~40-50%

However, GaAsP PMTs can be damaged by exposure to too much light.


Detectors - APDs



- Photons produce photoelectrons
- Photoelectrons are amplified by collisional ejection of electrons
- Small dynamic range




Hybrid Detectors




- GaAsP photocathode bombards an APD
- Very high gain
- Excellent for single photon counting
- Damaged by excessive light

Confocal vs. Widefield




Confocal




Widefield

Tissue culture cell with 60x / 1.4NA objective

Confocal vs. Widefield



Confocal




Widefield

20 μm rat intestine section recorded with 60x / 1.4NA objective

Doing more with your photons


Spectral Detection



Prism + Adjustable Slits: Leica SP microscopes

Doing more with your photons


Spectral Detection




Prism + PMT-array: Zeiss and Nikon

Doing more with your photons

Time-Correlated Single Photon Counting



The confocal microscope: Drawbacks



Scan excitation spot point-by-point to build up image

Problems:

Slow (~1 sec to acquire an image)

Low light efficiency (due to use of PMT as detector)

Solution:

Use multiple pinholes and a camera

A Solution: Spinning Disk Confocal







Image with many pinholes at once, so fast

Use CCD as detector, so much higher QE


The downside to the spinning disk



The downside to the spinning disk




The downside to the spinning disk



The downside to the spinning disk


Limited out of focus rejection




Pros/Cons of spinning disk

- Fast – multiple points are illuminated at once
- Photon efficient – high QE of CCD
- Gentler on live samples – usually lower laser power
- Fixed pinhole size (usually)
- Small field of view (usually)
- Crosstalk through adjacent pinholes limits sample thickness

Out-of-focus rejection



Examples



Drosophila S2 cell expressing
GFP-H2B and mCherry-tubulin
(Nico Stuurman and Ron Vale)



S. cerevisiae expressing a
mitochondrially targeted RFP,
Susanne Rafelski, Marshall lab

When to use confocal?

- Confocal is not a magic bullet
 - It is extremely wasteful of photons
 - Laser-scanning confocal is 100 – 200-fold less sensitive than widefield
 - Spinning-disk confocal is ~4-fold less sensitive than widefield

When to use confocal?

- Confocal is not a magic bullet
 - It is extremely wasteful of photons
 - High laser power generally result in more photobleaching and photodamage.
- For thin specimens, widefield epifluorescence is better – especially with deconvolution
- Confocal excels with thick, heavily stained specimens

When to use confocal?

- What is thick?
 - A good rule of thumb is $10 \times$ the depth of field of the objective
 - 100x / 1.4 NA: d.o.f. $0.66\mu\text{m}$
 - 20x / 0.75 NA: d.o.f. $2.3\mu\text{m}$
- Sample preparation is KEY for imaging thick specimens
 - Confocal does not fix scattering, refractive index mismatch, or everything else that can go wrong – it only removes out of focus light

TIRF

Total Internal Reflection Microscopy

Excite just the fluorophores at the coverslip.





Figure 1

Total Internal Reflection

Snell's Law: $n_1 \sin(\theta_1) = n_2 \sin(\theta_2)$

Aqueous sample
 $n_2 \approx 1.33\text{--}1.38$

Cover
Glass
 $n_1 = 1.518$




Total Internal Reflection

Snell's Law: $n_1 \sin(\theta_1) = n_2 \sin(\theta_2)$

Aqueous sample
 $n_2 \approx 1.33\text{--}1.38$

Cover
Glass
 $n_1=1.518$




Total Internal Reflection

Snell's Law: $n_1 \sin(\theta_1) = n_2 \sin(\theta_2)$

Aqueous sample
 $n_2 \approx 1.33\text{--}1.38$

Cover
Glass
 $n_1 = 1.518$




Total Internal Reflection

Snell's Law: $n_1 \sin(\theta_1) = n_2 \sin(\theta_2)$

Aqueous sample
 $n_2 \approx 1.33\text{--}1.38$

Cover
Glass
 $n_1=1.518$




Total Internal Reflection

Snell's Law: $n_1 \sin(\theta_1) = n_2 \sin(\theta_2)$

Aqueous sample
 $n_2 \approx 1.33\text{--}1.38$

Cover
Glass
 $n_1=1.518$




$$\sin(\theta_{\text{crit}}) = n_2 / n_1$$

Total Internal Reflection

Snell's Law: $n_1 \sin(\theta_1) = n_2 \sin(\theta_2)$

Aqueous sample
 $n_2 \approx 1.33\text{--}1.38$

Cover
Glass
 $n_1 = 1.518$





Total Internal Reflection

Decay length vs. angle


$$I \propto e^{-z/z_{TIRF}} = e^{-z \frac{4\pi}{\lambda} \sqrt{n_2^2 - n_1^2 \sin^2(\theta_1)}}$$

Typical TIRF depth
 $\approx 75\text{--}150\text{ nm}$




Two forms of TIRF microscopy

Prism coupled



- No excitation light in emission path
⇒ Very low background (if quartz slide)
- Needs separate, external beam path
⇒ Harder to align

Through the objective




- Easy to align
- Excitation light in emission path
⇒ Vulnerable to autofluorescence in the optics
- Requires very high NA

TIRF Objectives





Back focal plane
 $r \propto n \sin(\theta) \sim \text{"NA"}$



Typical NA 1.45–1.49

Extreme example: Olympus NA 1.65
Requires special high-index cover glass
and (volatile, toxic) immersion fluid

Epifluorescence vs. TIRF



Jaiswal et al 2002; cells loaded with FITC-dextran

Resources

<http://www.microscopyu.com>

<http://micro.magnet.fsu.edu>

James Pawley, Ed. “Handbook of Biological Confocal Microscopy, 3rd ed.”

http://www.hamamatsu.com/resources/pdf/etd/PMT_handbook_v3aE.pdf

Acknowledgements

Steve Ross, Mats Gustafsson