

(Image: T. Wittman, Scripps)

Monitoring molecular movement by microscopy

- General idea: optically mark one region of the cell, follow where it goes at later times.
- Implementation
 - Uncaging (old)
 - FRAP, FLIP: photobleaching
 - Photoactivation


FRAP

Fluorescence Recovery after Photo-bleaching




Need: to probe transport

Idea: bleach in one area,
watch recovery by transport from other areas

Fluorescence Recovery After Photobleaching (FRAP) with Green Fluorescent Protein



Measuring Cdc42 diffusion constant in yeast



Result: $df = (0.036 \pm 0.017) \mu\text{m}^2/\text{s}$


Marco et al. 2007 Cell 129:411-422

FLIP

Fluorescence Loss in Photo-bleaching

Need: probe connectivity

Idea: bleach in one compartment,
watch loss in connected compartments by exchange



Bleach one area repeatedly. Entire ER dims.
⇒ ER is contiguous


Photo-Activation

(Better?) FRAP/FLIP alternative

Some fluorophores can be activated by light

- Photo-uncagable dyes
- GFP-family proteins

Look for weak light
against dark background
Instead of slight dimming
of bright background



Activate a small area
Watch fluorescence *spread*

Photo-Activatable Proteins

Off-On


- PA-GFP, PS-CFP

Color change


- Kaede, KikGR, Eos,
- Dendra (activatable by blue)

Reversibly Switchable

- asCP, KFP (tetrameric)
- Dronpa






Dronpa – photoswitchable on and off



D R O N P A


Tracking actin flow with Dronpa



Hardware for photoactivation / photobleaching



- Need: controlled illumination of a region in the sample

Laser-scanning confocal microscope



Use mirrors to adjust angle at back focal plane – changes location of laser spot at sample


Laser-scanning confocal microscope




Laser-scanning confocal microscope

- Pros
 - can bleach / activate any shape ROI by scanning point by point
- Cons
 - Slow
 - Need a laser scanning confocal

Can control field of view with field




Better idea: Spatial Light Modulator





Can project any image onto sample

Digital Micromirror Device



Array of rotatable mirrors



Pros: Fast, can project any arbitrary image


Cons: Intensity is relatively low

FLIM

Fluorescence Lifetime Imaging

Measure the lifetime of the excited state


- Separate fluorophores with similar spectra
- Detect environmental parameters that affect lifetime (FRET, pH, ...)



FLIM


Measurement approaches

- Frequency domain
 - Modulated excitation
 - Lock-in detect emission phase




Time domain (pulsed exc.)

- Gated intensifier
 - Photon inefficient
- Time-correlated
 - single photon counting
 - Very efficient
 - \leq one photon per pulse \Rightarrow slow



FLIM



Examples





Hepatocyte membrane-stained with NBD, which has a hydrophobicity-dependent lifetime (TCSPC, 3 minutes for 300x300 pixels)

Environment-sensitive fluorophores



Many: pH, ions (Ca^{2+} , Mg^{2+} , Na^+ , K^+ , etc.), voltage, hydrophobicity,...



Fura-2








Fluorescein




MQAE

Environment-sensitive fluorophores




Imaging Cdc42 activation





Fuse sensor to GFP; I-SO/GFP ratio is proportional to degree bound

Monitoring Cdc42 activation



Fluorescence Resonance Energy Transfer




Sensitive to distances on the order of 10 nm, 20-fold below resolution limit

Distance dependence of FRET

$$E = \frac{1}{1 + (r^6/R_0^6)}$$

$$R_0^6 \propto k^2 n^{-4} QD J(l)$$

Overlap between donor
 emission and acceptor
 excitation
 Dipole-dipole interaction
 Donor quantum yield
 Refractive index




For CFP-YFP,
50% transfer at $R_0 = 4.9$ nm

Good FRET pairs


- CFP/YFP – use A206R mutants if dimerization is problematic
- GFP/mCherry or other FP pairs – not so well validated
- Fluorescein/Rhodamine
- Cy3/Cy5 or Rhodamine/Cy5
- Many other small molecule pairs

FRET Theory

- $k_2 = (\cos q_T - 3 \cos q_D \cos q_A)^2$
- For rapidly tumbling molecules, can average over all possible orientations to give $k_2 = 2/3$
- But rotational correlation time for GFP is ~ 16 ns; fluorescence lifetime is ~ 3 ns




Effects of FRET




- Donor lifetime shortened
- Acceptor emission depolarized
- Donor fluorescence quenched
- Acceptor fluorescence enhanced on donor excitation

Measuring FRET

- Donor lifetime shortened
- Can measure by fluorescence lifetime imaging, but requires specialized instrumentation




Measuring FRET




- Acceptor emission depolarized
- Can measure by fluorescence polarization microscopy

Measuring FRET




- Donor fluorescence quenched
- Acceptor fluorescence enhanced on donor excitation
- Can measure by donor recovery after acceptor photobleaching
 - Easy, but very sensitive to degree of photobleaching

Measuring FRET




- Donor fluorescence quenched
- Acceptor fluorescence enhanced on donor excitation
- Can measure by quantitative measurement of acceptor enhancement on donor excitation

A problem: crosstalk into FRET channel




Correct using
measurements from
CFP- and YFP- only cells


Crosstalk correction



Types of FRET experiments




Intramolecular



Intermolecular

Types of FRET experiments




Intramolecular

For intramolecular FRET, CFP and YFP are always present in a 1:1 ratio

Ratiometric imaging can be used as a rough measure of the amount of energy transfer

Types of FRET experiments





Intermolecular

For intermolecular FRET, the relative abundance of CFP and YFP is not controlled and can change over time.

Ratiometric imaging is no longer possible, and additional corrections are necessary.

Using FRET to monitor Rac activation



FCS

Fluorescence Correlation Spectroscopy


Small volume \rightarrow only a few molecules \rightarrow random fluctuations

- Study the *noise*
- Conclude about *random processes* at different *time scales*

Small excitation volume

Random molecular processes:

- Diffuse in and out
- Adopt different states
- Bind or react
- Photobleach
- ...



Study the auto-
(or cross-) correlation





Image Correlation Spectroscopy

Image FCS

Sense the random fluctuations of fluorescence within an *image* or *image sequence*

Much slower than point FCS, but get whole area
Can see where you are & deal with motion
Can do spatio-temporal analysis

Image series
Measure time and space
variations in intensity





Image Correlation Spectroscopy



Areas
analyzed

Flow
vector

○ = 10 min
diffusion
length



■ Stationary fraction
■ Flowing fraction
■ Diffusing fraction




Spatio-temporal
autocorrelation



Optical Tweezers

Mechanically manipulate the specimen with light

Why?



Measuring force and displacement of a single polymerase molecule



DNA recoil after stretching

Optical Tweezers


How does it work?

Photons carry momentum

Changing photon direction requires a force

More light refracted away from high intensity region

→ Force toward intensity maximum



Alternative way to look at it:



Field energy is less in higher index

System energy depends on particle position

$$dE/dx = \text{Force}$$

Holographic Optical Tweezers

Many traps at once
Independently movable
Made using a computer-controlled
spatial light modulator in a pupil plane



Further reading

www.microscopyu.com

micro.magnet.fsu.edu

Molecular Probes Handbook (probes.com)

James Pawley, Ed. "Handbook of Biological Confocal Microscopy, 3rd ed."

Acknowledgements

Mats Gustafsson