

A fluorescence microscopy image showing several elongated, curved structures, likely cells or tissue sections. The structures are stained with two different fluorescent dyes, resulting in a mix of green and red colors. The green signal is more prominent in some areas, while the red signal is more prominent in others. The overall image has a dark background, making the fluorescent structures stand out.

Live Cell Imaging

Key Elements for Successful Live Imaging

1. Choice of biological system
2. Validation
3. Reducing phototoxicity and photobleaching
4. Close preservation of normal environment

Choosing the Appropriate Organism to Answer your Question

1. Easy labeling of your target of interest
2. Biological process can be found at sufficient frequency
3. Biological process completes in a reasonable amount of time
4. The process can be sustained over the time it takes to image
5. Not imaging a process that occurs too deep into the tissue
6. The time resolution isn't beyond the imaging capabilities
7. Ways to validate your results?

Phototoxicity

Jablonski Energy Diagram




Figure 1


Biological Harmful Effects of Free Radicals

Because of its unpaired electron singlet O_2 will attack sites of high electron density -e.g. C=C, N atoms.

Lipid peroxidation


Protein modification

DNA modification




Yeast Viability After Imaging



10 msec exposures, 25 Z slices every 15 sec for 20 min



Using Meiotic Progression to Evaluate Phototoxicity



Tetrad Formation 3 fold decreased in imaged areas




Imaged Area



Not Imaged Area

Tetrads




>



Quantifying Phototoxicity in *C. elegans* Embryos

Image: 41 Z-slices every 2 min for 2 hours

Count number of nuclei




Tinevez et al. Meth. Enz.
506: 291 (2012)

Quantifying Phototoxicity in *C. elegans* Embryos

Image: 41 Z-slices every 2 min for 2 hours

Count number of nuclei



Tinevez et al. Meth. Enz.
506: 291 (2012)


Ways to Lessen Phototoxicity and Photobleaching

- Reduce signal (shorter exposures, lower intensity)
- Reduce frequency of exposures
- Remove oxygen and singlet oxygen
- Lower noise
- Improve detection

Dose Fractionation

- If cells can tolerate X amount of light:
- Can do 1000 exposures at $X/1000$
- 10 Z slices \times 100 time points
- 10 Z slices \times 4 colors \times 25 time points
- Etc.

Reduce frequency of exposures



Instead of 2x longer before bleaching, you can go 3x as long

1. Don't overwhelm endogenous cell antioxidant enzymes
2. Allow longer time for repair


Minimize Exposure

- Gate illumination by camera expose signal so that light is only on while camera is exposing
- Requires fast light source (laser or LED)

Minimize Exposure

- Gate illumination by camera expose signal so that light is only on while camera is exposing
- Requires fast light source (laser or LED)


Mechanical Shutter



Minimize Exposure

- Gate illumination by camera expose signal so that light is only on while camera is exposing
- Requires fast light source (laser or LED)

Fast Triggered Shutter



Removing Oxygen and Free Radicals

Oxygen scavengers



- Glucose oxidase / Catalase
- Protocatechuic acid/protocatechuate-3,4-dioxygenase
- Oxyrase (E. coli membrane particles)

Free radical scavengers

- Trolox – derivative of Vitamin E
- ascorbic acid – Vitamin C
- n-Propyl Gallate

GFP Bleaching in Live Cells


Incubation of cells with rutin 30 minutes prior to imaging increases photostability




Denoising Algorithms

- Kervrann – Patch based denoising (ND-SAFIR)
 - <http://serpico.rennes.inria.fr/doku.php?id=software:nd-safir:index>
- Luisier et al. wavelet denoising (SURE-LET)
 - <http://bigwww.epfl.ch/algorithms/denoise/>
- Also deconvolution

Patch-based denoising



SURE-LET denoising



Preserving a normal environment during imaging

Environmental Variables for Mammalian Cell Lines

Variable	Optimum Range	Comments
Temperature	28-37°C	Control with Specimen Chamber Heaters Inline Perfusion Heaters Objective Lens Heaters Environmental Control Boxes
Oxygenation	Variable	Perfuse or Change Media Regularly Use Large Chamber Volume
Humidity	97-100%	Closed chamber, humidified environmental chamber
pH	7.0 -7.7	Use Buffered Media, Perfuse or change media, no phenol red indicator
Osmolarity	260-320 mosM	Avoid evaporation, sealed chamber
Atmosphere	Air or 5-7% CO ₂	Use buffered media, closed chamber
Media buffer	Bicarbonate or Synthetic buffers	Beware of phototoxicity, closed & open chambers, atmosphere controlled chamber

Environmental Chamber







Figure 1


Bioptechs – Open/Closed Dish System



Delta T Dish



Heated Objective



Custom-made PDMS microfluidics

PDMS - Poly(dimethylsiloxane)


Relatively inexpensive

Easy to use Chemically inert/non-hazardous




Optically clear Flexible and fairly tough when cured

Easily bonded to itself or other materials


Permeable to air and liquids (but can be coated to prevent this).



CellASIC – Microfluidic Systems



CellASIC – Microfluidic Systems



Different Goals for Fixed vs. Live

Fixed specimens: Optimize the signal/noise

Live specimens: Optimize the signal/noise but without perturbing
your biological process

Typical ways to improve S/N in fixed samples

- Increase exposure time
- Amplify signal
- Brighter fluorophores/dyes
- Decrease photobleaching with antioxidants

>>>> Increase signal