



Superresolution Microscopy by Single Molecule Switching and Localization

Kurt Thorn

Nikon Imaging Center, UCSF



The Resolution Limit in Light Microscopy



- Image resolution improves with Numerical Aperture (NA)

$$NA = n \sin(\alpha)$$


where: α = light gathering angle
 n = refractive index of sample

The Resolution Limit in Light Microscopy



Improving The Resolution Limit in Light Microscopy

- $d = \lambda / 2n \sin \alpha$
- To improve resolution (smaller d):
- Use smaller λ – UV objectives
- Use larger n – Sapphire coverslips
- Use larger α – 4pi / I5M




Improving The Resolution Limit in Light Microscopy

- $d = \lambda / 2n \sin \alpha$
- What does this assume?
- Uniform illumination
- Continuous sample
- What if we break it?
- Structured illumination
- Localize single molecules (PALM/STORM/etc.)

How does imaging single molecules help us?


- Imagine taking a measurement
- Multiple times
- Adding measurements improves our ability to estimate the mean

$$SEM = \frac{\sigma}{\sqrt{N}}$$




How does imaging single molecules help us?

- Can measure the center of a distribution higher accuracy than any single measurement






How does imaging single molecules help us?


- Can measure the center of a distribution higher accuracy than it's width.
- Can determine the center of a single molecule to much higher accuracy than it's width.







Single-molecule localization precision



1 photon





$$d = \frac{\lambda}{2n \sin \alpha}$$




10 photons

100 photons

1000 photons




$$D \approx d / \sqrt{N}$$


Single molecules can be precisely localized

How can we apply this to imaging?


- Problem: in our image, we can't see single molecules – we see lots and lots of molecules and we can resolve them from each other.
- What if we could turn off all the molecules, and then turn them on a few at a time?

Super-resolution by localization


Fluorescence image



Raw images



STORM Image



2x real time

Stochastic Optical Reconstruction Microscopy = **STORM**

Also named as **PALM** (Betzig et al., Science, 2006) and **FPALM** (Hess et al., Biophys. J. 2006)




Reconstructed from 40,000 frames, 3,350,370 localization points

5 μm




B-SC-1 cell,
Microtubules stained with anti- β tubulin
Cy3 / Alexa 647 secondary antibody

Bates, Huang, Dempsey and Zhuang,
Science, 2007





500 nm

5 μm




5 μ m


A white horizontal scale bar representing 5 μ m.



3D Imaging: Localization in the Third Dimension





3D Imaging of the Microtubule Network



Huang, Wang, Bates and Zhuang,
Science, 2008

3D Imaging of the Microtubule Network



Small, isolated clusters



5 μm

Other 3D localization method


Astigmatic imaging



z (nm) -400 -200 0 200 400


Huang et al., Science 2008

Bi-plane imaging



Juette et al., Nat Methods 2008

Double-helical PSF



z (nm) -900 -500 0 600 1400

Pavani et al., PNAS 2009

Super-resolution by localization

- Combines two ideas:
 - Localization Microscopy
 - Single Molecule Switching

Single-Molecule Switching Approaches

- Lots of acronyms:

PALM

FPALM

STORM

dSTORM

GSDIM

PALMIRA

SMACM

PAINT

SPRAIPAINT




All use single-molecule
localization

Main differences are how you
switch the molecules on and
off


Single Molecule Switching using Fluorescent Proteins

PALM, FPALM

- Photoactivatable molecules: PA-GFP, PA-mCherry, PA-TagRFP



- Photoconvertible molecules: mEos2/3, PS-CFP2, PSmOrange




Parameters for evaluating dyes

- Number of photons emitted before bleaching / blinking



$$SEM = \frac{\sigma}{\sqrt{N}} \quad \Rightarrow \quad D = \frac{D_{diffraction}}{\sqrt{N_{photons}}}$$

- Contrast ratio – brightness of ON state relative to OFF state

Why Contrast Matters



- What if your OFF state isn't completely off?



- Can no longer distinguish single molecules above background


Fluorescent Proteins for Single Molecule Switching

Protein	λ_{ex}	λ_{em}	Contrast Ratio	Nphotons (median)	Notes
PS-CFP2	490	511	2000	260	Cyan-to-green
mGeos-M	503	514		387	Reversibly photoactivatable
PATagRFP	562	596	540	500?	Photoactivatable
PAmCherry	564	595	4000	348	Photoactivatable
tdEos	569	581		499	Green-to-red
mEos2	573	584		379	Green-to-red
mEos3.2	572	580		482	Green-to-red
PSmOrange	631	662	10700	337	Orange-to-far red


Full table at

<http://nic.ucsf.edu/dokuwiki/doku.php?id=storm:fps>


Multi-color and FPs




e.g. PA-GFP, mGeos-M




e.g. mEos2/3



e.g. PA-TagRFP, PA-mCherry



e.g. PS-CFP




e.g. PSmOrange

Single Molecule Switching Using Small Molecule Dyes

STORM, dSTORM, GSDIM

Singlet states

Triplet states




Single Molecule Switching Using Small Molecule Dyes

STORM, dSTORM, GSDIM

Singlet states

Triplet states




Single Molecule Switching Using Small Molecule Dyes

STORM, dSTORM, GSDIM

Singlet states

Triplet states

Other dark states




Dark states can be long-lived

Single Molecule Switching Using Small Molecule Dyes


STORM, dSTORM, GSDIM

- The idea:
- Drive the majority of molecules into a dark state
 - Use high laser power
 - Special buffers
- Wait for molecules to spontaneously return to ground state
 - Image them until the bleach or return to dark state


Photoswitching of red cyanine dyes



Multicolor Single Molecule Switching: Activation



Controlling the activation of Cy5



Cy3 / Alexa 647: Clathrin

Cy2 / Alexa 647: Microtubule

Crosstalk subtracted

Laser sequence


Cy3 — A647

Cy2 — A647




1 μ m

Multicolor Single Molecule Switching: Emission



- If you use high enough laser power, many molecules can be driven to the dark state
- Cy3B, Atto488, ...

Multicolor STORM/PALM



mEos2-tubulin



Alexa 647 anti- β tubulin

Drosophila S2 cells

2 μ m




Daichi Kamiyama



An alternative strategy: time slicing

- Use a multi-pass emission filter
- Alternate laser excitation so that every other frame is a different dye.



Multicolor imaging approaches

By emission wavelengths

- Simple fluorophores
- Low crosstalk
- Continuous imaging
- Multi-channel detection optics
- Needs nanometer scale image alignment

By activation wavelengths

- Dye-pairs
- Crosstalk from nonspecific activation
- Laser sequences
- Single channel detection
- Images naturally aligned

Summary: Labeling Strategies


- **Fluorescent proteins**
 - mEos is a good starting point
 - Multicolor may be challenging
- **Cy5 / A647 plus activator dyes**
 - Need to label your own antibodies
 - Crosstalk a problem
- **Dyes driven to dark state**
 - Need high laser power
- **Photoswitchable / photoactivatable dyes**
 - Not many options yet (Aberrior)

Effective resolution: Probe matters


Antibodies:
~ 10 nm

Fluorescent Proteins:
~ 3 nm


Small fluorophores:
~ 1 nm




Measured FWHM by antibody:	58 nm
Actual microtubule diameter:	25 nm
Measured FWHM by FP:	43 nm



~ 6000 photons



< 1000 photons



~ 6000 photons

Fluorescent protein vs. Antibody

Fluorescent protein fusion

- Live sample labeling
- High specificity
- High labeling efficiency
- Genetically encoded
- Lower S/N
- Multicolor imaging so far challenging

Antibody immunofluorescence

- Fixed sample
- Potential nonspecific labeling
- Lower labeling efficiency
- Labeling endogenous proteins
- High signal = high localization precision
- More versatile for multicolor imaging

Effective resolution: Density matters

Frames for image reconstruction:


200

500

1,000

5,000

40,000



Nyquist criteria



Point-to-point distance $\leq \frac{1}{2}$ Feature size

This labeling density limit of resolution applies to **all** fluorescence microscopy methods

Other issues

- Sample drift
 - Dedrifting by cross-correlation
 - Fiducial markers (gold particles, beads)
- Fixing
 - Background fluorescence
 - Sample preservation

Acknowledgements / Further Reading

- Bo Huang
- <http://nic.ucsf.edu/dokuwiki/doku.php?id=storm>
- Lots of papers