


# EMCCD Cameras

Rachit Mohindra  
Photometrics



## EMCCDs – The Need for Low Light Imaging


- Low excitation energy to reduce photo-bleaching/toxicity
  - Live cell fluorescence
  - Fast bleaching dyes
- Samples in ultra-low concentrations
  - $\sim 10^{-10}$  M in typical single molecule fluorescence applications
- Observe Processes with fast kinetics
  - Motility studies
  - Ion ratio imaging
- Essentially - Low Light and Speed

# CCD Image Acquisition

## Frame Transfer CCD




# EMCCD Cameras



**Back-illuminated, frame-transfer CCD (for >90% QE) or  
standard front-illuminated (65% QE)**

## EMCCD Cameras




## How Does this Help?

- If you are trying to measure signal that is lower than the read noise –  
You're not going to see anything
- If you multiply up the signal over the read-noise – Then you'll see something.
- But – the signal you're measuring isn't actually changing.
- You're effectively reducing the read noise


## Reduction of Read Noise

- Apply EM Gain factor of 200X
- Boost Signal 200X
- Equivalent to reducing the read noise by a factor of 200X
- If Read Noise is 50e-, after EM gain, it is equal to 0.25e-
- Allows extremely low-light detection

## Detect Signals Below the Read-Noise



Without EM Gain




With EM Gain



## Downfalls of EMCCDs

- Clock Induced Charge increases
  - Time independent (unlike dark current)
  - Increases background events
- Excess Noise Factor
  - Increases noise by a factor of  $\sim 1.4$
  - “Halves the Quantum Efficiency” – Not exactly
- Loss of Full Well Capacity
  - Inversely proportional to the amount of EM Gain applied
- Aging of EM Gain

## Excess Noise Factor





- Increases only Shot Noise and Dark Noise by a factor of  $\sim\sqrt{2}$  or  $\sim 1.4$
- Important at higher light levels, where images are dominated by Shot Noise

## Calibrate your EM Gain

- As EM Gain ages – setting 300X doesn't always provide 300X...
- Most Cameras provide you the ability to calibrate

## The Sweet Spot

- Maximum EM Gain is not the best setting
  - Maximizes Excess Noise
  - Maximizes Clock Induced Charge/Spurious noise
  - Drastically reduced Dynamic Range
  - Read Noise asymptotically reaches a minimum Read Noise
- Reduce the Read Noise to  $<1e^-$ 
  - No more EM Gain needed
  - Maximize the sensitivity, minimize the noise
  - Extend the life of the EM Gain



- With increasing EM Gain – ADU's increase, But Signal detected by the sensor remains constant
- After a certain point, the electron count will remain the same

## EMCCD's – A Summary

- Low Light Imaging that requires Speed
- Effectively reduces your EM Gain
- Not a perfect technology
- Find the Sweet Spot – Don't need more EM than this

# Thank You!

For more information please visit:-

[www.photometrics.com](http://www.photometrics.com)