



A fluorescence microscopy image showing several cells. Some cells exhibit a bright red signal, while others show a bright green signal. The background is dark, and there are some smaller, dimmer spots of red and green light scattered throughout the field of view.

FRET and Biosensors


Kurt Thorn
Nikon Imaging Center

Image: Thomas Huckaba


Fluorescence Resonance Energy Transfer




Cameleons: FRET-based Ca^{2+} sensors



Cameleons: FRET-based Ca^{2+} sensors





Calcium transients in *C. elegans* pharynx



Kerr et al. 2000. *Neuron* 26, p. 583-594

Using FRET to monitor Rac activation



Good FRET pairs


- CFP/YFP – use A206R mutants if dimerization is problematic
- GFP/mCherry, YFP/mCherry, mTFP/mKO, many other FP pairs – not so well validated
- Fluorescein/Rhodamine
- Cy3/Cy5 or Rhodamine/Cy5
- Many other small molecule pairs

Distance dependence of FRET

$$E = \frac{1}{1 + (r^6/R_0^6)}$$




$$R_0^6 \propto \kappa^2 n^{-4} Q_D J(\lambda)$$

Overlap between donor
 emission and acceptor
 excitation
 Donor quantum yield
 Refractive index
 fluorophores



For CFP-YFP,
 50% transfer at $R_0 = 4.9$ nm

Transition dipole of GFP



Angular dependence of FRET

$$R_*^6 \propto n^{-4} Q_D J(\lambda)$$

$$E = \frac{1}{1 + (r^6 / R_*^6 \kappa^2)}$$


κ^2 depends on the relative orientations of the donor and acceptor excitation dipoles.

κ^2 ranges between 0 and 4 and is 0 for whenever the donor and acceptor dipoles are perpendicular to one another.

For rapidly-rotating dyes
 $\kappa^2 = 2/3$


FRET Theory

- $\kappa^2 = (\cos \theta_T - 3 \cos \theta_D \cos \theta_A)^2$
- For rapidly tumbling molecules, can average over all possible orientations to give $\kappa^2 = 2/3$
- But rotational correlation time for GFP is ~ 16 ns; fluorescence lifetime is ~ 3 ns




Effects of FRET

- Donor lifetime shortened
- Acceptor emission depolarized
- Donor fluorescence quenched
- Acceptor fluorescence enhanced on donor excitation




Measuring FRET

- Donor lifetime shortened
- Can measure by fluorescence lifetime imaging, but requires specialized instrumentation




Measuring FRET

- Acceptor emission depolarized
- Can measure by fluorescence polarization microscopy




Measuring FRET

- Donor fluorescence quenched
- Acceptor fluorescence enhanced on donor excitation
- Can measure by donor recovery after acceptor photobleaching
 - Easy, but very sensitive to degree of photobleaching




Measuring FRET


- Donor fluorescence quenched
- Acceptor fluorescence enhanced on donor excitation
- Can measure by quantitative measurement of acceptor enhancement on donor excitation



Types of FRET experiments





Intramolecular



Intermolecular

Types of FRET experiments




Intramolecular

For intramolecular FRET, CFP and YFP are always present in a 1:1 ratio

Ratiometric imaging can be used as a rough measure of the amount of energy transfer

Types of FRET experiments




Intermolecular

For intermolecular FRET, the relative abundance of CFP and YFP is not controlled and can change over time.


Ratiometric imaging is no longer possible, and additional corrections are necessary.

Data Acquisition


Three things to measure:



Donor
Intensity



FRET
Intensity



Acceptor
Intensity


Data Acquisition

- Maximize signal-to-noise: use high NA objective, sensitive, low-noise camera, high-transmission filters
- Minimize shifts between wavelengths
 - Fluor or apochromatic objective
 - Multipass dichroic with external excitation and emission filters

Image preprocessing


- Background subtraction
- Register images by maximizing correlation with FRET image

Data Acquisition




Acquire sequential images of FRET, YFP, CFP, and DIC

A problem: crosstalk into FRET channel



Correct using measurements from CFP- and YFP- only cells

A problem: crosstalk into FRET channel



For strains with only CFP and YFP,
 $\text{FRET}_C = 0$

Fit


$$\text{FRET}_C = \text{FRET}_m - \alpha \text{CFP} - \beta \text{YFP} - \gamma$$

Typical values:


$$\alpha \sim 0.9$$

$$\beta \sim 0.4$$

Crosstalk correction




Calculating FRET efficiency



Traditionally:

$$E = 1 - \frac{F_D (\text{Donor+Acceptor})}{F_D (\text{Donor alone})}$$


Calculating FRET efficiency




$$E = 1 - \frac{FRET_C \cdot G + F_D}{F_D}$$

G corrects for detection efficiencies of CFP and YFP


$$G = Q_D \Phi_D / Q_A \Phi_A$$



FRET



FRETc




Efficiency

One final issue: Autofluorescence

We correct for autofluorescence in the FRET channel by inclusion of γ


But we also need to correct for autofluorescence in the donor channel



$$E = \frac{\text{FRET}_C \cdot G}{\text{FRET}_C \cdot G + F_D}$$

Correct donor autofluorescence by subtracting median donor fluorescence of untagged cells

Data analysis procedure



Preprocessing:

Background subtraction


Image alignment by maximizing the correlation of donor and acceptor with the FRET image.

Typical shifts are <2 pixels

Photobleaching


- Some dyes photobleach quite easily (prime offenders: fluorescein, YFP)
- Correction procedures are available but are non-trivial
- Photobleaching can lead to peculiar artifacts

Spatial variation of efficiency



Illumination Uniformity


HBO 103



0.96 0.98 1 1.02 1.04




HBO 103 + LLG



0.96 0.98 1 1.02 1.04



Lambda LS + LLG



0.96 0.98 1 1.02 1.04





FRET Conclusions


- Using FRET as a qualitative reporter is relatively straightforward.
- Quantitative FRET is challenging and requires correction of a large number of potential artifacts.
- Trying to use FRET to infer distances *in vivo* is probably best avoided.
- Choice of fluorescent proteins for FRET is likely to be idiosyncratic and system dependent.

Single domain sensors – GCaMP2


Intensity-based sensor



Single domain sensors – GCaMP2





Ca^{2+} bound




Ca^{2+} free

Single domain sensors - roGFP


Ratiometric sensor



Imaging glutathione redox potential *in vivo*




Coupling glutaredoxin-1 to roGFP makes it specifically sensitive to glutathione redox potential and accelerates its response



Bimolecular fluorescence complementation






BiFC (aka split GFP)




Has speed and reversibility issues,
so most useful as a screening tool.

Environment-sensitive fluorophores

Cdc42-binding domain of WASP as sensor for active Cdc42





Imaging Cdc42 activation



Fuse sensor to GFP; I-SO/GFP ratio is proportional to degree bound

Monitoring Cdc42 activation



Additional reading

- Lakowicz, “Principles of Fluorescence Spectroscopy”, Chapters 13-15
- Gordon et al. 1998, Biophys. J. **74** p. 2702-2713
- Berney and Danuser 2003, Biophys. J. **84** p.3992-4010
- Zal and Gascoigne 2004, Biophys. J. **86** p 3923-3939
- FRET code is at: <https://github.com/kthorn/fretproc>