




(Image: T. Wittman, Scripps)

The Light Microscope

- Four centuries of history
- Vibrant current development
- One of the most widely used research tools





Major Imaging Functions of the Microscope

- Magnify
- Resolve features
- Generate Contrast
- Capture and Display Images

An Upright Epifluorescence Microscope






Figure 1

Electromagnetic Waves




Waves vs. Photons vs. Rays

- Quantum wave-particle duality
- Rays: photon trajectories
- Rays: propagation direction of waves




Rays are perpendicular to wavefronts



Light travels more slowly in matter

The speed ratio is the ***Index of Refraction, n***

$$v = c/n$$




Refractive Index Examples

• Vacuum	1
• Air	1.0003
• Water	1.333
• Cytoplasm	1.35–1.38 ?
• Glycerol	1.475 (anhydrous)
• Immersion oil	1.515
• Fused silica	1.46
• Optical glasses	1.5–1.9
• Diamond	2.417

Depends on wavelength and temperature

Refraction by an Interface



⇒ Snell's law:

$$n_1 \sin(\theta_1) = n_2 \sin(\theta_2)$$

Mirror law:


$$\theta_r = \theta_1$$

Which Direction?


n_1

$n_2 > n_1$

Refraction goes
towards the normal
in the *higher-index* medium




Lenses work by refraction




Ray Tracing Rules of Thumb

(for thin ideal lenses)


Parallel rays converge
at the focal plane




Rays that cross in the focal plane
end up parallel



Rays through the lens center are unaffected



Imaging







The lens law:

$$\frac{1}{L_1} + \frac{1}{L_2} = \frac{1}{f}$$

Magnification:

$$M = \frac{d_2}{d_1} = \frac{L_2}{L_1}$$


Real and virtual images



The same lens law applies: Negative lenses have negative f
 Virtual objects or images have negative values of $L1$ or $L2$


Finite vs. Infinite Conjugate Imaging

- Finite conjugate imaging (older objectives)




Infinite conjugate imaging (modern objectives). Image at infinity

⇒ Need a **tube lens**



Magnification: $M = \frac{f_1}{f_0}$


Back focal plane




Rays that leave the object with the same
angle

meet in the objective's *back focal plane*


The Compound Microscope




The Compound Microscope




The Compound Microscope




The Compound Microscope



The Compound Microscope



The Compound Microscope



Eyepieces (Oculars)


Aberration-Free 10x Eyepiece With Diopter Adjustment




Features

Magnification (10x typical)
“High eye point” (exit pupil high enough to allow eyeglasses)
Diopter adjust (at least *one* must have this)
Reticle or fitting for one
Eye cups

Trans-illumination Microscope




Köhler Illumination




- Each light source point produces a parallel beam of light at the sample
- Uniform light intensity at the sample even if the light source is "ugly" (e.g. a filament)

Critical Illumination



- The source is imaged onto the sample
- Usable only if the light source is perfectly uniform

Conjugate Planes in A Research Microscope




How view the pupil planes?

Two ways:

- “Eyepiece telescope”
- “Bertrand lens”


By far the most important part:
the Objective Lens



Each major manufacturer sells 20-30 different ***categories*** of objectives.
What are the important distinctions?

Working Distance

Objective Working and Parfocal Distance



In general, high NA lenses have short working distances

However, extra-long working distance objectives do exist


Some examples:

10x/0.3 WD = 15.2mm

20x/0.75 WD = 1.0mm


100x/1.4 WD = 0.13mm

The focal length of a lens depends on the refractive index...




$$f \propto 1/(n-1)$$

... and the refractive index depends on the wavelength (“dispersion”)




⇒ **Chromatic aberration**

- Different colors get focused to different planes
- Not good...




Dispersion vs. refractive index of different glass types




Achromatic Lenses

- Use a weak negative flint glass element to compensate the dispersion of a positive crown glass element




Achromats and Apochromats




Correction classes of objectives



Achromat
t
(cheap)



Fluor
“semi-apo”
(good
correction,
high UV
transmission)



Apochromat
(best
correction)

**Correction for other (i.e. monochromatic) aberrations
also improves in the same order**—————→

Curvature of Field


Focal plane
Focal surface

Tube lens

objective


sample

Focal
surface



Plan objectives

- Corrected for field curvature
- More complex design
- Needed for most photomicrography



- **Plan-Apochromats** have the highest performance (and highest complexity and price)

Putting one brand of objectives onto another brand of microscope?

Usually a bad idea:

- May not even fit
 - May get different magnification than is printed on the objective
 - Incompatible ways of correcting lateral chromatic aberration (LCA)
- ⇒ mixing brands can produce severe LCA



Tube lens focal length


Nikon	200
Leica	200
Olympus	180
Zeiss	165

LCA correction:

<u>In objective</u>	<u>In tube lens</u>
Nikon	Leica
Olympus	Zeiss


Lateral chromatic aberration

(= **LCA, lateral color, chromatic difference of magnification**)
= **Different magnification for different colors**




Interference


In phase



constructive interference




Opposite phase





destructive interference



Diffraction by a periodic structure (grating)



Diffraction by a periodic structure (grating)




In phase if:

$$d \sin(\theta) = m \lambda$$


for some integer m

Diffraction by an aperture

drawn as waves




Light spreads to new angles





Larger aperture
 \Leftrightarrow
weaker diffraction

Diffraction by an aperture

drawn as
rays



The pure, “far-field”
diffraction pattern
is formed at ∞ distance...





...or can be formed
at a finite distance
by a lens...


...as happens in a microscope

The Airy Pattern


= the far-field diffraction pattern from a round aperture




"Airy disk" diameter
 $d = 2.44 \lambda f/d$
(for small angles d/f)




Aperture and Resolution




Diffraction spot
on image plane
= ***Point Spread
Function***




Aperture and Resolution




Diffraction spot
on image plane
= ***Point Spread
Function***




Aperture and Resolution




Diffraction spot
on image plane
= ***Point Spread
Function***



Aperture and Resolution



Diffraction spot
on image plane
(resolution)




- Image resolution improves with aperture size ~~Aperture size~~ Numerical Aperture (NA)


$$NA = n \sin(\alpha)$$

where: α = light gathering angle
 n = refractive index of sample

Numerical Aperture



4X / 0.20 NA
 $\alpha = 11.5^\circ$




100X / 0.95 NA
 $\alpha = 71.8^\circ$

Numerical Aperture

Compare:


Numerical Aperture:

$$NA = n \sin(\alpha)$$




Snell's law:

$$n_1 \sin(\theta_1) = n_2 \sin(\theta_2)$$



- $n \sin(\theta)$ doesn't change at horizontal interfaces
 - $\sin(\text{anything}) \leq 1$
- ⇒ NA cannot exceed the *lowest* n between the sample and the objective lens




Numerical Aperture

Compare:


Numerical Aperture:

$$NA = n \sin(\alpha)$$




Snell's law:

$$n_1 \sin(\theta_1) = n_2 \sin(\theta_2)$$




- $n \sin(\theta)$ doesn't change at horizontal interfaces
- $\sin(\text{anything}) \leq 1$

- ⇒ NA cannot exceed the *lowest* n between the sample and the objective lens
- ⇒ $NA > 1$ requires **fluid immersion**



Immersion Objectives



NA can approach
the index of the
immersion fluid

Oil immersion:

$n \approx 1.515$

max NA \approx **1.4** (1.45–1.49
for TIRF)

Glycerol immersion:

$n \approx 1.45$ (85%)


max NA \approx **1.35** (Leica)

Water immersion:

$n \approx 1.33$

max NA \approx **1.2**

Objective Types



Basic properties

- Magnification
- Numerical Aperture (NA)
- Infinite or finite conjugate
- Cover slip thickness if any
- Immersion fluid if any

Correction class

- Achromat
- Fluor
- Apochromat

Field flatness

- Plan or not

Phase rings for phase contrast

- Positive or negative
- Diameter of ring (number)

Special Properties

- Strain free for Polarization or DIC

Features

- Correction collar for spherical aberration
- Iris
- Spring-loaded front end
- Lockable front end

Further reading

www.microscopyu.com

micro.magnet.fsu.edu

Douglas B. Murphy “Fundamentals of Light Microscopy and Electronic Imaging”

James Pawley, Ed. “Handbook of Biological Confocal Microscopy, 3rd ed.”

Acknowledgements

Ron Vale / Mats Gustafsson